

ANAESTHETICS

ANAESTHETICS



\$33.55
each

Anaesthetics Cartridge Holder

Holds 50 1.8ml anaesthetics cartridges.
Size: 89mm (Diameter) x 76mm (H)
Cleaning instructions:
Use only soap and water on acrylic surface. Alcohol or ammonia-based disinfectants may cause cracks to appear.

Code	Colour
993830	Clear



\$39.15
each




Needle Station

Holds any dental needle sheath.
Made of stainless steel with silicone non-slip ring base and fully autoclavable up to 135°C.
Size: 50mm (Diameter) x 38mm (H)

Code	Colour
993831	Metallic

NuMedical



\$10.95
/box

100pcs/box

Needle Sheath Prop

A disposable needle recapping device that can be used to safely uncap and recap needles.
Size: 70mm (L) x 70mm (H)

Code
993832

ANAESTHETICS



\$29.95
each

For 2.2mL anesthetic injection



Syringe, Aspirating, 2.2mL, Point Type

Used for injection of anaesthetics.

Code	Colour
996794	Metallic

ANAESTHETICS



Nu2medical

ARTICULATING

ARTICULATING



\$36.95 /box

300pcs/box



Bausch Articulating Paper, W/Dispenser

200U(0.008"). Easily obtain accurate pressure distribution within seconds. High spots are immediately visible.



Made in Germany

Code	Size	Colour
993091	BK-01	Blue
993092	BK-02	Red



\$33.95 /box

300pcs/box



Bausch Articulating Paper, Refillbox

200U(0.008"). Easily obtain accurate pressure distribution within seconds. High spots are immediately visible.



Made in Germany

Code	Size	Colour
993093	BK-1001	Blue
993094	BK-1002	Red



\$17.50 /box

200pcs/box



Bausch Articulating Paper, Pre-cut Strips

40U(0.0016"). Easily obtain accurate pressure distribution within seconds. High spots are immediately visible.



Made in Germany

Code	Size	Colour
993095	BK-61	Blue
993096	BK-62	Red



\$29.95 /box

200pcs/box



Bausch Articulating Paper, Booklet

40U(0.0016"). Easily obtain accurate pressure distribution within seconds. High spots are immediately visible.



Made in Germany

Code	Size	Colour
993097	BK-09	Blue
993098	BK-10	Red



\$27.95 /box

50pcs/box



Bausch Articulating Paper, Horseshoe

200U(0.008"). Easily obtain accurate pressure distribution within seconds. High spots are immediately visible.



Made in Germany

Code	Size	Colour
993088	BK-03	Blue
993089	BK-04	Red



\$36.95 /box

300pcs/box

Bausch Articulating Paper, Booklet

200U(0.008"). Easily obtain accurate pressure distribution within seconds. High spots are immediately visible.



Made in Germany

Code	Size	Colour
993090	BK-05	Blue

Nu2medical

ARTICULATING



\$17.95
/box



Articulating Paper, Horseshoe (100 microns)

Easily obtain accurate pressure distribution within seconds. High spots are immediately visible.

Code	Quantity	Colour
993080	20pcs/book, 6 books/box	Red
993085	12pcs/book, 6 books/box	Blue/Pink



\$6.55
/box



Articulating Paper, Booklet (100 microns)

Easily obtain accurate pressure distribution within seconds. High spots are immediately visible.

Code	Quantity	Colour
993083	20pcs/book, 10 books/box	Blue
993084	12pcs/book, 12 books/box	Blue/Red



\$6.55
/box

20pcs/book
10 books/box



Articulating Paper, Booklet (30 microns)

Easily obtain accurate pressure distribution within seconds. High spots are immediately visible.

Code	Colour
993081	Red



\$6.55
/box


20pcs/book
10 books/box



Articulating Paper, Booklet (30 microns)

Easily obtain accurate pressure distribution within seconds. High spots are immediately visible.

Code	Colour
993082	Blue



\$9.95
each

Articulating Paper Forceps, Straight

Used to hold articulating paper when checking occlusion of dental restorations and teeth.

Code
996750



\$9.95
each

Articulating Paper Forceps, Curved

Used to hold articulating paper when checking occlusion of dental restorations and teeth.

Code
996751

ARTICULATING

NuMedical