

# COMPOSITE & CEMENT

COMPOSITE & CEMENT

\$16.95  
/box

6pcs/box



Scan for details  
and cautions

**POLISHING KIT FOR COMPOSITE,  
INTRAORAL USE**

Code	Size
991224	RA0306E

\$57.95  
/box

12pcs/box



Scan for details  
and cautions

**POLISHING KIT FOR COMPOSITE,  
INTRAORAL USE**

Code	Size
991223	RA0312D

\$63.95  
/box

9pcs/box



Scan for details  
and cautions

**POLISHING KIT FOR PORCELAIN,  
EXTRAORAL USE**

Code	Size
991222	HP0209D

\$57.95  
/box

12pcs/box



Scan for details  
and cautions

**POLISHING KIT FOR PORCELAIN,  
INTRAORAL USE**

Code	Size
991221	RA0212D

\$157.95  
/box

9pcs/box




Scan for details  
and cautions

**POLISHING KIT FOR ZIRCONIA,  
EXTRAORAL USE**

Code	Size
991220	HP0109D

\$109.95  
/box

12pcs/box




Scan for details  
and cautions

**POLISHING KIT FOR ZIRCONIA,  
INTRAORAL USE**

Code	Size
991219	RA0112D

NuMedical

# COMPOSITE & CEMENT



\$157.95  
/box

10pcs/box



Scan for details  
and cautions

## POLISHING KIT FOR GLASS - CERAMICS, EXTRAORAL USE

Code	Size
991218	HP1210D



\$109.95  
/box

12pcs/box



Scan for details  
and cautions

## POLISHING KIT FOR PRESSIBLE PORCELAIN, INTRAORAL USE

Code	Size
991217	RA1212D



\$34.75  
/box

## Polishing Kit

Polishing Kit with 50x Super Fine, 50x Fine, 50x Medium and 50x Coarse polishing discs, and 4 mandrels.

Code
993144



\$15.65  
/box

50pcs/box

## POLISHING DISC REFILLS

Code	Size	Colour
993145	Super Fine	White
993146	Fine	Green
993147	Medium	Pink
993148	Coarse	Black



\$239.00  
each

## Amalgamator YG-100, Amalgam and Capsule

Amalgamator YG-100 ensures a precise and consistent mixing of amalgam and capsules (Note: capsules that are 15mm in length are compatible). The amalgamator deactivates automatically once the front cover is opened. Mixing time: 1 - 99 secs

Code
993107

COMPOSITE & CEMENT

NuMedical

# COMPOSITE & CEMENT

COMPOSITE & CEMENT

Nu2medical



\$6.55  
each

Composites  
not included.

## Syringe Organiser

Holds up to 20 composites.  
Size: 220mm (L) x 120mm (W) x 40mm (H)

Code	Type
994179	Stand Set



\$3.95  
each

## Syringe Organiser

Code	Type
994177	Insert

Buy 3 get 1 free



\$53.95  
/bottle

3ml/bottle

## DMP Single Bonding

One component light curing adhesive material for the bonding of enamel and dentin in aesthetic restorations. It combines primer and adhesives into a solution to be applied after etching enamel and dentin simultaneously (total etch wet bonding technique). Solvent-free material and compatible with all composites. Included: 3ml Bottle x 1.

Code
991216

Buy 3 get 1 free



\$12.95  
each

## DMP Etching Gel - 3ml

An acid etchant for enamel and dentin. Its phosphoric acid content is 37%. Etching gel increases the surface area of the tooth by making it rougher thus resulting in better adhesion for the bonding agent. For ease of dispensing the syringe tip is bendable so that any angle required can be obtained. Included: 3ml Syringe x1 + tips.

Code
991215

Buy 3 get 1 free



\$52.35  
each

## NANOCERAM Bright Flow Composite

NANOCERAM Bright Flow is a light curing flowable version of the Nanoceram-Bright composite using the same nanohybrid technology to achieve optimum results. It is characterized by its superior aesthetics through its chameleon effect and high gloss polishability. This highly thixotropic product guarantees easy handling and precise application. Included: 2g Syringe x2

Size	Code
A2	991210

Buy 3 get 1 free



\$52.35  
each

## NANOCERAM Bright Composite

NANOCERAM Bright is a light cure nano-hybrid composite composed of a special blend of micro and nano fillers which produces a superior restorative. The highest aesthetic results have been achieved through the development of a sophisticated method for the integration of the nanoparticles in the resin matrix. These nanoparticles, of average particle size of 20 nm, assist in a quick polish after placement as well as ultralow wear. Included: 4g Syringe x1

Size	Code
A2	991169
A3	991170

# COMPOSITE & CEMENT

COMPOSITE & CEMENT



\$3.95  
/bottle

100pcs/bottle

## Micro Applicators, Regular, 100mm L

Used to apply cavity liners, bonding agents, etchants, fluoride, hemostatic and disclosing solutions. Tips are made from non-absorbent, lint-free fibre. Single-use. Available in light and dark colours which are useful as an indicator for different applications. Tip Size: 2.5mm. Length: 100mm

Code	Colour
990801LB	L. Blue



\$3.95  
/bottle

100pcs/bottle

## Micro Applicators, Regular, 100mm L

Used to apply cavity liners, bonding agents, etchants, fluoride, hemostatic and disclosing solutions. Tips are made from non-absorbent, lint-free fibre. Single-use. Available in light and dark colours which are useful as an indicator for different applications. Tip Size: 2.5mm. Length: 100mm

Code	Colour
990801DB	D. Blue



\$3.95  
/bottle

100pcs/bottle

## Micro Applicators, Fine, 100mm L

Used to apply cavity liners, bonding agents, etchants, fluoride, hemostatic and disclosing solutions. Tips are made from non-absorbent, lint-free fibre. Single-use. Available in light and dark colours which are useful as an indicator for different applications. Tip Size: 2.0mm. Length: 100mm

Code	Colour
990803LG	L. Green



\$3.95  
/bottle

100pcs/bottle

## Micro Applicators, Fine, 100mm L

Used to apply cavity liners, bonding agents, etchants, fluoride, hemostatic and disclosing solutions. Tips are made from non-absorbent, lint-free fibre. Single-use. Available in light and dark colours which are useful as an indicator for different applications. Tip Size: 2.0mm. Length: 100mm

Code	Colour
990803DG	D. Green



\$3.95  
/bottle

100pcs/bottle

## Micro Applicators, Ultra Fine, 100mm L

Used to apply cavity liners, bonding agents, etchants, fluoride, hemostatic and disclosing solutions. Tips are made from non-absorbent, lint-free fibre. Single-use. Available in light and dark colours which are useful as an indicator for different applications. Tip Size: 1.5mm. Length: 100mm

Code	Colour
990805LP	L. Purple



\$3.95  
/bottle

100pcs/bottle

## Micro Applicators, Ultra Fine, 100mm L

Used to apply cavity liners, bonding agents, etchants, fluoride, hemostatic and disclosing solutions. Tips are made from non-absorbent, lint-free fibre. Single-use. Available in light and dark colours which are useful as an indicator for different applications. Tip Size: 1.5mm. Length: 100mm

Code	Colour
990805DP	D. Purple



# COMPOSITE & CEMENT

COMPOSITE & CEMENT



\$3.95  
/bottle

100pcs/bottle

## Micro Applicators, Cylinder, 100mm L

Used to apply cavity liners, bonding agents, etchants, fluoride, hemostatic and disclosing solutions. Tips are made from non-absorbent, lint-free fibre.

Single-use. Available in light and dark colours which are useful as an indicator for different applications.  
Tip Size: 1.2mm. Length: 100mm

Code	Colour
990807	Black




\$18.65  
/box

60pcs/box

## Sticky Applicators

Applicators designed for picking up small objects, colour-transfer and tissue-marking. Locating denture high spots just got easier. Unique, flexible dual-bend handle design to reach tight areas. Single-use.

Code	Colour
990843	White




\$27.45  
each

To release applicators,  
turn the handle  
clockwise as per  
indicated on  
the dispenser.

## Micro Applicator Dispenser (One-Way)

Capacity: 400pcs  
Size: 141mm (L) x 100mm (W) x 90mm (H)  
Cleaning instructions:  
Use only soap and water on acrylic surface.  
Alcohol or ammonia-based disinfectants may cause cracks to appear.

Code
990808




\$37.65  
each

## Micro Applicator Dispenser (Two-Way)

Capacity: 400pcs  
Size: 125mm (L) x 155mm (W) x 90mm (H)  
Cleaning instructions:  
Use only soap and water on acrylic surface.  
Alcohol or ammonia-based disinfectants may cause cracks to appear.

Code
990809



## Mixing Wells

Compatible with all dental materials. Easy to pick up. Single-use.

Type	Quantity	Price	Code	Colour
2 Wells	200pcs/bag	\$6.25/bag	990790	White
4 Wells	200pcs/bag	\$7.65/bag	990791	White

NuMedical

# COMPOSITE & CEMENT

## COMPOSITE & CEMENT

NuMedical



\$0.85  
each

### Universal Brush Handle

Made of durable and lightweight plastic for repeated autoclaving. Suits all composite brushes and brush tips. Autoclavable up to 135°C. Single-use.

Code	Colour
990838	White
990839	Blue
990840	Green
990841	Red
990844	Yellow



\$6.55  
/bag

100pcs/bag

### Brush Tips

Tapered tip for smooth, even application of all types of liquid. These brush tips fits our Universal Brush Handle and most brush handles. Single-use.

Code  
990842



100mm

\$7.65  
/bottle

100pcs/bottle

### Brush Applicators

Durable brush which can be bended to any angle. Easy to snap apart for lingual or interproximal use. Comes in assorted colours of blue, yellow, red and green. Single-use. Size: 100mm (L)

Code  
990816 Assorted



\$93.50  
each

### Tri-Spectra Shade Guide With Led

Efficient shade-taking with light modes; daylight mode to determine initial shade, room light and ambient light modes to verify shade.

Code  
993087



\$94.55  
/box

### Shade Guide

For comparing teeth shades before and after bleaching only. R-20 Bleaching Shade Guide.

Code  
993086



\$14.95  
/box

50pcs/box

### Monojet Syringe - Curved

Designed for surgical site and post-surgical home irrigation. Curved tip to reach tight areas. Can be used with silicone and rubberbased impression materials. Single-use. Size: 12cc(12ml).

Code  
991097

# COMPOSITE & CEMENT

COMPOSITE & CEMENT






\$5.50  
/bag

100pcs/bag

## Etch Tips

Designed for placing etchants, sealants, cements and other flowable materials. Tips feature a blunt end and are pre-bent to the appropriate angle for convenience. Single-use.

Size	Code	Colour
19g (Composites, cements)	990866	Black
18g (Cavity liners, bases, composites)	990867	Pink
20g (Composites, cements)	990868	Yellow
22g (Sealants, etchants, cements)	990869	Grey
25g (Etchants)	990870	Blue



## Mixing Pads - Paper + Poly

100pcs/book, 3 books/box

Poly-coated for superior performance. Single-use.



Size	Price	Code
45mm (L) x 35 mm (W)	\$4.90/box	990851
80mm (L) x 70mm (W)	\$8.10/box	990852
140mm (L) x 80mm (W)	\$10.95/box	990854

**Buy 4 get 1 free**

## Non-Skid Mixing Pads - Paper + Poly

100pcs/book, 3 books/box


Poly-coated and non-skid poly-laminated to prevent slipping. Single-use.



Size	Price	Code
80mm (L) x 70mm (W)	\$10.95/box	990853
140mm (L) x 80mm (W)	\$12.50/box	990855

**Buy 4 get 1 free**


NuMedical



## Mixing Pads - Plastic


100pcs/book, 3 books/box

These pads are made of plastic-coated special paper, which can be easily detached from the pad. Single-use.



Size	Price	Code
80mm (L) x 70mm (W)	\$8.75/box	990856
140mm (L) x 80mm (W)	\$12.50/box	990858


**Buy 4 get 1 free**



## Non-Skid Mixing Pads - Plastic

100pcs/book, 3 books/box

Plastic-coated and non-skid poly-laminated to prevent slipping. Single-use.



Size	Price	Code
80mm (L) x 70mm (W)	\$12.50/box	990857
140mm (L) x 80mm (W)	\$13.95/box	990859

**Buy 4 get 1 free**



**\$24.30 each**

## Carpules Dispenser Gun

Ergonomic design, Suitable for most capsules, Fast loading, For all viscosity materials, Autoclavable up to 135°C.

Code	Colour
990836	Blue
990837	White



# COMPOSITE & CEMENT


COMPOSITE & CEMENT



## Multi-Purpose Dappen Dish

Suitable for prophy, etchant and bonding. Single-use. Size: 21mm (H) x 20mm (Diameter)

Code	Colour
990793	Assorted



## Glass Dappen Dish

Size: 30mm (H) x 22mm (Diameter)  
Compatible with all dental materials.  
Assorted colours of blue yellow white and green. *Cleaning instructions: Wipe with alcohol-based disinfectant solution.*


Code	Colour
990795	Assorted



## Universal Materials Well


Suitable for prophy etchant and bonding. Autoclavable up to 135°C. Size: 31mm (H) x 18mm (Diameter)

Code	Colour
990887	Blue
990888	Green
990889	Orange
990890	Purple
990891	White




## UV Multipurpose Holder

Size: 326mm (L) x 195mm (W) x 50mm (H)  
*Cleaning instructions: Use only soap and water on acrylic surface. Alcohol or ammonia-based disinfectants may cause cracks to appear.*



Code	Colour
996973	Orange



## UV Cover

For protection of light-curing materials from ultraviolet light.  
*Cleaning instructions: Use only soap and water on acrylic surface. Alcohol or ammonia-based disinfectants may cause cracks to appear.*


Size Options:  
Small: 100mm (L) x 75mm (W) x 30mm (H)  
Large: 185mm (L) x 115mm (W) x 30mm (H)

Size	Price	Code
Small	\$10.30 each	990876
Large	\$12.50 each	990877

# COMPOSITE & CEMENT

## COMPOSITE & CEMENT

**Nu**medical



**\$65.95 each**


Accessories not included.

### Composites and Microbrush Holder, Type 1

12 syringe holes + 4 Compartments + 1 Microbrush Holder.  
Size: 222mm (L) x 170mm (W) x 96mm (H).  
*Cleaning instructions:*  
Use only soap and water on acrylic surface. Alcohol or ammonia-based disinfectants may cause cracks to appear.



Code	Colour
996970	Clear



**\$98.95 each**


Accessories not included.

### Composites and Microbrush Holder, Type 2

20 syringe holes + 4 Compartments + 2 Microbrush Holders.  
Size: 334mm (L) x 170mm (W) x 96mm (H).  
*Cleaning instructions:*  
Use only soap and water on acrylic surface. Alcohol or ammonia-based disinfectants may cause cracks to appear.



Code	Colour
996972	Clear



**\$89.95 each**

Accessories not included.

### Composites and Microbrush Holder, Type 3

20 syringe holes + 8 Compartments + 1 Microbrush Holder.  
Size: 328mm (L) x 170mm (W) x 96mm (H).  
*Cleaning instructions:*  
Use only soap and water on acrylic surface. Alcohol or ammonia-based disinfectants may cause cracks to appear.



Code	Colour
996971	Clear



**\$17.55 each**


Accessories not included.

### Composite Holder - 12 Cross Shaped

12 cross-shaped syringe holes. Size: 118mm (L) x 118mm (W) x 68mm (H).  
*Cleaning instructions:*  
Use only soap and water on acrylic surface. Alcohol or ammonia-based disinfectants may cause cracks to appear.



Code
990862



**\$28.75 each**

Accessories not included.

### Composite Holder - 24 Bell Shaped

24 bell-shaped holes. Size: 220mm (L) x 118mm (W) x 68mm (H).  
*Cleaning instructions:*  
Use only soap and water on acrylic surface. Alcohol or ammonia-based disinfectants may cause cracks to appear.



Code
990863



**\$17.55 each**

Accessories not included.

### Composite Holder - 12 Round

12 round syringe holes. Size: 118mm (L) x 118mm (W) x 68mm (H).  
*Cleaning instructions:*  
Use only soap and water on acrylic surface. Alcohol or ammonia-based disinfectants may cause cracks to appear.



Code
990850

# COMPOSITE & CEMENT

COMPOSITE & CEMENT


\$28.75  
each



Accessories  
not included.


## Composite Holder - 24

24 round syringe holes.  
Size: 220mm (L) x 118mm (W) x 68mm (H).  
*Cleaning instructions:*  
Use only soap and water on acrylic surface.  
Alcohol or ammonia-based disinfectants may cause cracks to appear.



Code  
990835

\$25.25  
each



Accessories  
not included.

## Multi-Purpose Acrylic Holder, Type 1

6 syringe holes + 4 Compartments.  
Size: 100mm (L) x 145mm (W) x 85mm (H).  
*Cleaning instructions:*  
Use only soap and water on acrylic surface.  
Alcohol or ammonia-based disinfectants may cause cracks to appear.



Code  
990830

Colour  
Clear

\$35.15  
each



Accessories  
not included.

## Multi-Purpose Acrylic Holder, Type 2


14 syringe holes + 8 Compartments.  
Size: 200mm (L) x 145mm (W) x 85mm (H).  
*Cleaning instructions:*  
Use only soap and water on acrylic surface.  
Alcohol or ammonia-based disinfectants may cause cracks to appear.



Code  
990831

Colour  
Clear

\$46.15  
each



Accessories  
not included.

## Multi-Purpose Acrylic Holder, Type 3


30 syringe holes + 12 Compartments.  
Size: 294mm (L) x 172mm (W) x 95mm (H).  
*Cleaning instructions:*  
Use only soap and water on acrylic surface.  
Alcohol or ammonia-based disinfectants may cause cracks to appear.



Code  
990834

Colour  
Clear

\$15.35  
each



Accessories  
not included.

## Composite Holder - 11


11 round syringe holes.  
Size: 100mm (L) x 130mm (W) x 110mm (H).  
*Cleaning instructions:*  
Use only soap and water on acrylic surface.  
Alcohol or ammonia-based disinfectants may cause cracks to appear.



Code  
990832

Colour  
Clear

\$23.80  
each



Accessories  
not included.

## Composite Holder - 23

23 round syringe holes.  
Size: 200mm (L) x 130mm (W) x 110mm (H).  
*Cleaning instructions:*  
Use only soap and water on acrylic surface.  
Alcohol or ammonia-based disinfectants may cause cracks to appear.



Code  
990833

Colour  
Clear



## COMPOSITE & CEMENT

COMPOSITE & CEMENT



Nu<sup>o</sup>medical