



Highly Durable Files
Uniquely Designed
Strong Cutting Force
Rebound Memory

Exclusively Distributed By





\$29.95
/pack

UDG M3 Large Taper Gold NiTi Rotary Files

Tapered design with gold coating to shape canals.

(*ProTaper Gold™ Users. ProTaper Gold™ is a registered trademark of Dentsply.)

6pcs/pack

Size	LENGTH			
	21mm	25mm	31mm	
S1	995096	995101	995157	
S2	995097	995102	995158	
F1	995098	995103	995159	
F2	995099	995104	995160	
F3	995100	995105	995161	
Assorted (SX, S1, S2, F1, F2, F3)	995151	995152	995162	

ENDODONTICS

NuMedical



UDG M3 Path NiTi Rotary

For optimal shaping of the root canal.

(*Tulsa Dental Path File™ Users.

Tulsa Dental Path File™ is a registered trademark of Dentsply.)

\$29.95
/pack

6pcs/pack

Size	LENGTH	
	21mm	25mm
Assorted (2pcs each size)	995106	995107

ENDODONTICS

ENDODONTICS



\$29.95
/pack

6pcs/pack



UDG M3 Pro Gold NiTi Rotary Files

Refined variable tapered path files with gold coating by UDG. (*ProTaper Gold™ Users. ProTaper Gold™ is a registered trademark of Dentsply.)

SIZE							
Length	#17/12	#12/02	#16/02	#18/05	#25/06	#35/04	#17/12 - #35/04
19mm	995067						
21mm	-	995068	995069	995070	995071	995072	995066
25mm	-	995075	995076	995077	995078	995079	995073

SIZE						
Length	#20/04	#20/06	#30/04	#30/06	#35/06	#40/06
25mm	995369	995370	995371	995372	995373	995374

Nu2medical



**SPECIALLY
DESIGNED FOR
BABY TEETH!**

\$29.95
/pack

6pcs/pack

UDG M3 Immature NiTi Rotary Files

Designed with a shorter length for endo treatments on baby teeth. (*Researched and developed by UDG.)

LENGTH	
Size	16mm
Assorted (2pcs each size)	995080



\$31.95
/bag

10pcs/bag

Endo Training Block

Size: 10mm (L) x 10mm (W) x 30mm (H)

Code	Colour
995108	Clear



UDG M3-L Platinum One File System

Continuous innovation of the Rotary One-File System by UDG.
(*WaveOne™ Gold and VDW Reciproc™ Users. WaveOne™ Gold is a registered trademark of Dentsply. VDW Reciproc™ is a registered trademark of VDW.)

\$29.95
/pack

6pcs/pack

Length	SIZE			
	#20/07	#25/06	#35/06	#20/07 - #35/06
21mm	995083	995084	995085	995081
25mm	995086	995087	995088	995082

M3-LPLATINUM

One File System

ISO Colour Coding

M3-L Platinum instruments are ISO colour-coded for easy identification

Depth Marking

M3-L Platinum instruments have radiograph visible depth markings.
Eg: 18, 19, 20 and 22mm.

Silicone Stopper

The stopper is ISO colour-coded.
Yellow: 0.20mm
Red: 0.25mm
Green: 0.35mm

Shaft

M3-L Platinum instruments have a short shaft of 11mm to enable better access to molars compared to many other instruments which have a shaft of 13mm or longer.

Patented flat-sided S-shaped cross section for better and more efficient debris removal.

Working Length

Working Length Depth Markings at:


- 21mm
- 25mm
- 31mm
- 18mm, 19mm and 20mm
- 18mm, 19mm, 20mm and 22mm
- 22mm and 24mm

NuMedical

ENDODONTICS

ENDODONTICS

ENDODONTICS



UDG Reamer Hand Files, Stainless Steel

Suitable for safe root canal preparation with reaming motion.



\$5.95
/pack

6pcs/pack

SIZE									
Length	#6	#8	#10	#15	#20	#25	#30	#35	#40
21mm	995248	995141	995142	995143	995144	995145	995249	995250	995251
25mm	995260	995146	995147	995148	995149	995150	995261	995262	995263
31mm	995272	995273	995274	995275	995276	995277	995278	995279	995280

SIZE								
Length	#45	#50	#55	#60	#70	#80	#15-#40	#45-#80
21mm	995253	995254	995255	995256	995257	995258	995252	995259
25mm	995265	995266	995267	995268	995269	995270	995264	995271
31mm	995282	995283	995284	995285	995286	995287	995281	995288

NuMedical




UDG Hand File Holder

Ideal for hand files with better treatment of posterior teeth and patients with small mouth opening. Medical class stainless steel use. Autoclavable up to 135 °C.

\$59.95
each

Code
995307



UDG M3-C Hand Files, Stainless Steel

Suitable for severely curved or calcified root canals.

\$6.95
/pack

6pcs/pack

SIZE						
Length	#6	#8	#10	#12.5	#15	#6-#15
19mm	995111	995112	995113	995114	995115	995168
21mm	995116	995117	995118	995119	995120	
25mm	995163	995164	995165	995166	995167	

At 18, 19, 20 & 22mm

Inactive Pilot Tip (conducts the instrument safely along the canal)




\$5.95
/pack

UDG K Hand Files, Stainless

Suitable for all preparation techniques with filing motion.

6pcs/pack



SIZE									
Length	#6	#8	#10	#15	#20	#25	#30	#35	#40
21mm	995169	995121	995122	995123	995124	995125	995170	995171	995172
25mm	995181	995126	995127	995128	995129	995130	995182	995183	995184
31mm	995193	995194	995195	995196	995197	995198	995199	995200	995201

SIZE								
Length	#45	#50	#55	#60	#70	#80	#15-#40	#45-#80
21mm	995174	995175	995176	995177	995178	995179	995173	995180
25mm	995186	995187	995188	995189	995190	995191	995185	995192
31mm	995203	995204	995205	995206	995207	995208	995202	995209

ENDODONTICS



\$5.95
/pack

6pcs/pack

UDG M3-H Hand Files, Stainless

Suitable for all preparation techniques with filing motion.
(ProTaper Gold™ Users. ProTaper Gold™ is a registered trademark of Dentsply.)

Nu
medical

SIZE								
Length	#8	#10	#15	#20	#25	#30	#35	#40
21mm	995131	995132	995133	995134	995135	995210	995211	995212
25mm	995136	995137	995138	995139	995140	995221	995222	995223
31mm	995232	995233	995234	995235	995236	995237	995238	995239

SIZE								
Length	#45	#50	#55	#60	#70	#80	#15-#40	#45-#80
21mm	995214	995215	995216	995217	995218	995219	995213	995220
25mm	995225	995226	995227	995228	995229	995230	995224	995231
31mm	995241	995242	995243	995244	995245	995246	995240	995247

NIC[®]

U⁺ NiTi Files - ProTaper Gold Users

For opening and preparation of the root canal.

\$39.95
/pack

6pcs/pack

Size	LENGTH			
	19mm	21mm	25mm	
SX	995560			
S1		995473	995481	
S2		995474	995482	
F1		995475	995483	
F2		995476	995484	
F3		995477	995485	
F4		995479	995487	
F5		995480	995488	
SX - F3		995478	995486	 S1 S2 F1 F2 F3 SX



\$7.65
/pack

6pcs/pack

NiTi K Files

Suitable for all preparation techniques with filing motion.

SIZE								
Length	#8	#10	#15	#20	#25	#30	#35	#40
21mm	993670	993671	993672	993673	993674	993675	993676	993677
25mm	993653	993654	993655	993656	993657	993658	993659	993660

SIZE								
Length	#45	#50	#55	#60	#70	#80	#15 - #40	#45 - #80
21mm	993678	993679	993680	993681	993682	993683	993684	993685
25mm	993661	993662	993663	993664	993665	993666	993667	993668



\$5.95
/pack

6pcs/pack

K Files, Stainless Steel

Suitable for all preparation techniques with filing motion.

	SIZE								
Length	#6	#8	#10	#15	#20	#25	#30	#35	#40
21mm	993601	993602	993603	993604	993605	993606	993607	993608	993609
25mm	993584	993585	993586	993587	993588	993589	993590	993591	993592

SIZE								
Length	#45	#50	#55	#60	#70	#80	#15 - #40	#45 - #80
21mm	993610	993611	993612	993613	993614	993615	993616	993617
25mm	993593	993594	993595	993596	993597	993598	993599	993600



\$7.65
/pack

6pcs/pack

NiTi H Files

For efficient enlargement of root canals in filing motion.

SIZE								
Length	#8	#10	#15	#20	#25	#30	#35	#40
21mm	993619	993620	993621	993622	993623	993624	993625	993626
25mm	993636	993637	993638	993639	993640	993641	993642	993643

SIZE								
Length	#45	#50	#55	#60	#70	#80	#15 - #40	#45 - #80
21mm	993627	993628	993629	993630	993631	993632	993633	993634
25mm	993644	993645	993646	993647	993648	993649	993650	993651



\$5.95
/pack

6pcs/pack

H Files, Stainless Steel

For efficient enlargement of root canals in filing motion.

	SIZE								
Length	#6	#8	#10	#15	#20	#25	#30	#35	#40
21mm	993550	993551	993552	993553	993554	993555	993556	993557	993558
25mm	993567	993568	993569	993570	993571	993572	993573	993574	993575

SIZE								
Length	#45	#50	#55	#60	#70	#80	#15 - #40	#45 - #80
21mm	993559	993560	993561	993562	993563	993564	993565	993566
25mm	993576	993577	993578	993579	993580	993581	993582	993583

ENDODONTICS



H Rotary Files, Stainless Steel

For efficient enlargement of root canals in filing motion.

\$10.65
/pack

6pcs/pack

Length	SIZE						
	#15	#20	#25	#30	#35	#40	#15 - #40
21mm	995496	995497	995498	995499	995500	995501	995502
25mm	995503	995504	995505	995506	995507	995508	995509



Paste Carriers, Stainless Steel

For filling root canals with paste, cement or sealer.

\$11.40
/pack

4pcs/pack

Length	SIZE				
	#25	#30	#35	#40	#25 - #40
21mm	995456	995457	995458	995459	995460
25mm	995461	995462	995463	995464	995465



Peeso Reamers, Stainless Steel

For widening and straightening the coronary part of the root canal.

\$13.95
/pack

6pcs/pack

Length	SIZE						
	#1	#2	#3	#4	#5	#6	#1 - #6
28mm	993686	993687	993688	993689	993690	993691	993698
32mm	993692	993693	993694	993695	993696	993697	993699



Gates Drills, Stainless Steel

For enlargement of the coronal or the middle third of the root canal.

\$13.95
/pack

6pcs/pack

Length	SIZE						
	#1/50	#2/70	#3/90	#4/110	#5/130	#6/150	#1 - #6
32mm	995489	995490	995491	995492	995493	995494	995495

ENDODONTICS



NiTi Spreaders

Suitable for lateral condensation in curved canals.

\$7.65
/pack

6pcs/pack

Length	SIZE						
	#15	#20	#25	#30	#35	#40	#15 - #40
21mm	995428	995429	995430	995431	995432	995433	995434
25mm	995435	995436	995437	995438	995439	995440	995441



Spreaders, Stainless Steel

Suitable for lateral condensation in curved canals.

\$5.95
/pack

6pcs/pack

Length	SIZE						
	#15	#20	#25	#30	#35	#40	#15 - #40
21mm	995442	995443	995444	995445	995446	995447	995448
25mm	995449	995450	995451	995452	995453	995454	995455



NiTi Pluggers

Suitable for vertical condensation of gutta-percha cones.

\$7.65
/pack

6pcs/pack

Length	SIZE						
	#15	#20	#25	#30	#35	#40	#15 - #40
21mm	995400	995401	995402	995403	995404	995405	995406
25mm	995407	995408	995409	995410	995411	995412	995413



Pluggers, Stainless Steel


Suitable for vertical condensation of gutta-percha cones.

5.95
/pack

6pcs/pack

Length	SIZE						
	#15	#20	#25	#30	#35	#40	#15 - #40
21mm	995414	995415	995416	995417	995418	995419	995420
25mm	995421	995422	995423	995424	995425	995426	995427

ENDODONTICS



\$5.70
/pack

10pcs/pack

Barbed Broaches

For easy and safe removal of pulp tissue.

SIZE								
Length	#0/XXXXF	#1/XXXF	#2/XXF	#3/XF	#4/F	#5/M	#6/C	All Sizes
21mm	995466	995467	995468	995469	995470	995471	995472	995553



\$29.95
/pack

6pcs/pack

M3 Large Tapered NiTi Hand Files

Suitable for all preparation techniques with filling motion.
(Protaper Hand File Users. Protaper is a registered trademark of Dentsply.)

SIZE						
Length	S1	S2	F1	F2	F3	Assorted
21mm	995289	995290	995291	995292	995293	995294
25mm	995295	995296	995297	995298	995299	995300
31mm	995301	995302	995303	995304	995305	995306



\$81.95
/pack

Endo Pluggers, Set of 4

Each set comes with 1/1 0.5mm yellow, 1/2 0.5mm/0.6mm red, 2/3 0.6mm/0.8mm blue and 2/4 0.6mm/0.1mm green.

Code
996845

ENDODONTICS

Nu
medical

Gutta Percha Points

- Saves time with easy filling for faster and easier root canal obturation.
- Reaches the APX and completely seals the length of root canal.
- Length marking available for points of 18, 19, 20, 22 and 24mm.

		Size								
	Pack	Price	#15	#20	#25	#30	#35	#40	#45	#15 - #40
0.02	120pcs	\$6.95	995513	995514	995515	995516	995517	995518	995519	995525
0.04	60pcs	\$6.95	995527	995528	995529	995530	995531	-	-	995539
0.06	60pcs	\$7.95	-	995541	995542	995543	-	-	-	995552




		Size	
	Pack	Price	
XF	100pcs	\$9.85	995008
FF	100pcs	\$9.85	995009
M	100pcs	\$9.85	995011

Paper Points

- Quick and easy depth measurements.
- Completely dries the full length of the root canal.
- Length marking available for points of 18, 19, 20, 22 and 24mm.

	Size								
	Pack	Price	#15	#20	#25	#30	#35	#40	#15 - #40
0.02	200pcs	\$6.95	995561	995562	995563	995564	995565	995566	995573
0.04	100pcs	\$7.25	995575	995576	995577	995578	995579	995580	-
0.06	100pcs	\$7.95	-	995589	995590	995591	-	-	995000



		Size	
	Pack	Price	
F	200pcs	\$8.45	995003
M	200pcs	\$8.45	995004
C	200pcs	\$8.45	995005

ENDODONTICS



\$97.60
each

Endo Multi Holder

Ideal for cleaning and autoclaving endo instruments. Holds 16 files (all lengths). Autoclavable up to 135°. Size: 65mm (L) x 28mm (W) x 60mm (H). Note: Files not included.

Code
993001



\$5.95
/pack

FREE! UDG Endo File Sample Box, Please Ask Your Regional Sales For Getting Free Endo Kit.

Ideal for customer trying our endo file products. It includes 6 pieces of C UDG Endo Hand File, 6 pieces of K UDG Endo Hand File, 6 pieces of UDG M3 Pro Gold Rotary Endo Files and 3 pieces of Endo Training Block.

Code
995308



ONE YEAR
WARRANTY



View Details

Easydo Activator

The activator provides a solution for the defects and deficiencies of the current clinically commonly used ultrasonic root canal cleaning devices and instruments that are insufficient in the apex cleaning efficiency, easily cause root canal wall perforation and the risk of cross-infection. The Activator with a portable device improve the treatment effect. Safe, Portability and Easy of Using.

\$359.00
/pack

Code
993700



View Details

Activator Needle Tips

Designed to be compatible with 993700. The tip is soft and good toughness and can be bent. It can reduce and avoid the risk of breakage and not damage the root canal that has been dredged. The vortex is formed in the high-frequency oscillation, causing the liquid to produce tens of thousands of tiny air pressure bubbles. These bubbles form an energy burst and explode, impact the inner wall of root canal and root tip, and quickly peel off dirt, impurities, and debris. At the same time, paste can be delivered and the broken needle can be taken out.


\$129.00
/pack

60pcs/pack

Code
993701

ENDODONTICS

NuMedical



\$12.95 each


Instruments not included.

Endodontic Case - Large D

Ideal for keeping up to 30 H Files and reamers. Sludge on files can be easily removed by the file cleaner. The gauze can be replaced. Autoclavable up to 135°. Size: 107mm (L) x 57mm (W) x 52mm (H)



Code	Colour
993043	Ivory



\$12.95 each

Instruments not included.

Bur Case - Large E

Ideal for autoclave sterilisation. Suitable for FG, RA and HP bur. Autoclavable up to 135°. Size: 107mm (L) x 57mm (W) x 52mm (H)



Code	Colour
993042	Ivory



\$12.95 each

Endodontic Case - Large A

Ideal for holding H Files and reamers with a total of 60 holes. Autoclavable up to 135°. Size: 107mm (L) x 57mm (W) x 52mm (H). Note: Files not included.



Code	Colour
993046	Ivory




\$12.95 each

Endodontic Case - Large C

Ideal for keeping H Files, reamers, gutta percha and paper points. Can be sterilised at the same time. Autoclavable up to 135°. Size: 107mm (L) x 57mm (W) x 52mm (H)



Code	Colour
993044	Ivory



\$12.95 each

Instruments not included.

Endodontic Case - Large B

Ideal for keeping colour-coded gutta percha points and paper points. Autoclavable up to 135°. Size: 107mm (L) x 57mm (W) x 52mm (H)



Code	Colour
993045	Ivory



\$6.95 each

Endodontic Case - Small


Ideal for storing cotton rolls, cotton balls or small instruments. Can be sterilised. Autoclavable up to 135°. Size: 57mm (L) x 57mm (W) x 57mm (H)



See our cotton rolls Pg 58

Code	Colour
993040	Ivory

ENDODONTICS



\$9.95
each


Instruments
not included.



Endodontic Case - Small F

Ideal for keeping H-Files and reamers.
Autoclavable up to 135°.
Size: 57mm (L) x 57mm (W) x 57mm (H)

Code	Colour
993065	Ivory



\$9.95
each


Instruments
not included.



Endodontic Case - Small G

Sludge on files can be easily removed by the file cleaner. The gauze can be replaced.
Autoclavable up to 135°.
Size: 57mm (L) x 57mm (W) x 57mm (H)

Code	Colour
993066	Ivory



\$9.95
each

Instruments
not included.



Endodontic Case - Small H

Ideal for keeping gutta-percha and paper points classified by colour.
Autoclavable up to 135°.
Size: 57mm (L) x 57mm (W) x 57mm (H)

Code	Colour
993067	Ivory




\$8.75
/bag

100pcs/bag

Endo Stops

These are easy-to-use silicone, round measuring stops that help establish the canal depth and measurement of each file during root canal therapies. Autoclavable up to 135°C.

Code	Colour
993012	Red
993013	White
993014	Yellow
993074	Dark Blue
993075	Blue
993076	Black
993077	Purple
993078	Green




\$15.25
each

Endo Ruler - Stainless Steel

Made from 304 food-grade stainless steel. Ruler scale provides markings up to 35mm in 1mm increments. Ruler can be used by right and left-handed users. Autoclavable up to 135°C.
Size: 38mm (L) x 11mm (W)

Code	Colour
993068	White



\$7.65
each


Endo Ruler

Ruler scale provides marking up to 35mm in 1mm increments. Ruler can be used by right and left-handed users. Autoclavable up to 135°C.
Size: 45mm (L) x 15mm (W)

Code	Colour
993009	Blue

ENDODONTICS

Nu
medical



\$7.65 each

NuEndo Holder Type 1

This NuEndo Holder holds 11 Engine and 5 Hand files with measurement gauge, sponge insert and stopper. Right or left hand use. Autoclavable up to 135°C. Contains Antimicrobial Protection. Each set comes with 1 NuEndo Holder and 12 sponge inserts. Size: 138mm (L) x 10mm (W) x 55mm (H)

Code	Colour	Code	Colour
995031	White	995039	Red
995032	Blue	995040	Navy Blue
995033	Green	995041	Copper
995034	Beige	995042	Yellow
995035	Gray	995043	Navy Green
995036	Plum	995044	Orange
995037	Teal	995045	Purple
995038	Dark Blue	995046	Pink



\$7.65 each

NuEndo Holder Type 2

This NuEndo Holder holds 16 Hand files or Reamers with measurement gauge, sponge insert and stopper. Right or left hand use. Autoclavable up to 135°C. Contains Antimicrobial Protection. Each set comes with 1 NuEndo Holder and 12 sponge inserts. Size: 138mm (L) x 10mm (W) x 55mm (H)

Code	Colour	Code	Colour
995047	White	995055	Red
995048	Blue	995056	Navy Blue
995049	Green	995057	Copper
995050	Beige	995058	Yellow
995051	Gray	995059	Navy Green
995052	Plum	995060	Orange
995053	Teal	995061	Purple
995054	Dark Blue	995062	Pink



\$7.65 each




NuEndo Holder Type 3

This NuHolder Holds 8 files with measurement gauge and stopper. 2 Sterilisation counters on locking bar tracks the number of times instruments have been sterilised. Autoclavable up to 135 degree. Size: 75mm(L) x 55mm(W).

Code	Colour
996394	White
996395	Blue
996396	Red
996397	Green
996398	Purple

\$9.95 each



Cutta-Percha Tip Trimming and Cutting

Measure and cut the Gutta Percha Points. Autoclavable up to 135 degree.

Code
995154

ENDODONTICS



\$5.65
each

Cutta Percha Points Measuring Ruler

Measure the size of Gutta Percha Point. Autoclavable up to 135 degree.

Code
995153



\$6.55
/bag

48pcs/bag

NuEndo Sponge Inserts

Sponge inserts for NuEndo Ring and NuEndo Holders. Disposable.
Size: 45mm (L) x 22mm (W) x 15mm (H)

Code
995063

Colour
Black



\$6.55
/bag

48pcs/bag

Sponge Inserts for NuEndo Stand

Sponge inserts for NuEndo Stand. Disposable.
Size: 57mm (L) x 57mm (W) x 5mm (H)

Code
993004

Colour
Black



\$4.35
/set

NuEndo Stand

NuEndo Stand contains Antimicrobial Protection properties and holds Endo instruments in an upright position. Thicker foam prevents solution from leaking with rounded corners for easy cleaning. Autoclavable up to 135°. Each set comes with 1 NuEndo Stand and 13 foam refills.
Size: 45mm (L) x 45mm (W) x 35mm (H)

Code
993007

Colour
Black



\$6.55
/set

NuEndo Ring

NuEndo Ring is manufactured with Antimicrobial Protection with built-in millimeter gauge. Suitable for both right and left hand use. Autoclavable up to 135°C. Each set comes with 1 NuEndo Ring and 13 foam refills.
Size: 39mm (L) x 40mm (W) x 32mm (H)

Code	Colour
995015	White
995016	Blue
995017	Green
995018	Beige
995019	Gray
995020	Plum
995021	Teal
995022	Dark Blue


Code	Colour
995023	Red
995024	Navy Blue
995025	Copper
995026	Yellow
995027	Navy Green
995028	Orange
995029	Purple
995030	Pink

ENDODONTICS

NuMedical

ENDODONTICS

ENDODONTICS



\$9.95 /set

Endo Stand

Made from high-heat resistant plastic. Used to clean, sterilise and stage instruments. Each set includes 1 Endo Stand and 13 sponge inserts. Autoclavable up to 135°C. Size: 35mm (H) x 55mm (Diameter)

Code	Colour
993003W	White
993003B	Blue
993003G	Green
993003R	Pink
993003P	Purple



\$9.95 /set

Triangular Endo Ring

Made from high-heat resistant plastic. Used to clean sterilise and stage instruments. Autoclavable up to 135°C. *Cushion included.


Code	Colour
993002	Blue




\$17.80 each

Endo Measuring Block - Plastic

Made from high-heat resistant plastic. Easy-to-use and accurate. Autoclavable up to 135°. Size: 100mm (L) x 24mm (W) x 32mm (H)



Code	Colour
993008	Ivory



\$12.25 /box

50pcs/box

Endo Cushions - Round

Disposable endo cushions. Compatible with Endo Stand (#993003). Size: 65mm (Diameter)

Code	Colour
993006	Blue



Luer-Lock Syringe

For surgical site irrigation. Single-use.

Size	Pack	Price	Code
3ml	100pcs/box	\$8.55/box	993056
6ml	100pcs/box	\$8.55/box	993057
12ml	100pcs/box	\$12.90/box	993058



\$28.55 /box

50pcs/box



Notched-tip

Sterile Irrigation Kit

Luer-lock syringes with notched-tip needles (30g) for irrigating root canals.

Size	Code
3ml	990897
6ml	990898
12ml	990899

ENDODONTICS

ENDODONTICS

Nu²medical



\$21.95
/bag

50pcs/bag

Root Canal Needle

Pre-bent needle tip to access and clean root canals.

Code	Colour
995065	Clear



\$16.55
/bag

25pcs/bag

HVE Luer Adaptor

For suction of root canal irrigation, sodium hypochlorite, acid etching and adhesive agent, etc.
Fit for different types of syringe needles for various treatments.




Code	Colour
995064	Clear



\$21.95
/box

20pcs/box



Double-Vented

Irrigating Needle Tips


Double-vented irrigating needle tips for excellent irrigation. The needle tips are bendable.

Size	Code	Colour
22mm, 30g	995109	Yellow
26mm, 30g	995110	Blue



\$23.40
/box

100pcs/box



Single-Vented

Irrigating Needle Tips - Closed End

Bendable needle tip at 1°. Used to irrigate subgingival surgical sites and the sulcus. Closed-end with side vent for use in apical proximity. Universal luer-style design. Single-use.

Size	Code	Colour
23g	993052	Blue
27g	993053	Yellow
30g	993054	Purple



\$14.70
/bag

50pcs/bag

Mixing Tips, Brown/White+Tip, 1:1

Disposable mixing tips ideal for mixing dental cold-whitening cements, resins, curing cements, self-curing adhesive materials, crown and bridge materials, composite inlays, etc.

Code	Colour
991087	Brown




\$14.70
/bag

50pcs/bag

Mixing Tips, Brown/White+Tip, 1:1

Disposable mixing tips ideal for mixing dental cold-whitening cements, resins, curing cements, self-curing adhesive materials, crown and bridge materials, composite inlays, etc.

Code	Colour
991091	Brown



\$359.00
/pack

ONE YEAR
WARRANTY

Code
992929

NuMedical R7 Apex

1.3cm Ultra-Thin Body Design with Type-C Measurement Cable/Power Connector

Host shape with stylish ultra-thin design, beautiful, simple, atmospheric.
New type-C measurement cable/power cable interface, convenient for daily use, reject a variety of interfaces.

3.5" TFT LCD Screen

Accurate visual experience
320x480 resolution, color saturation and vividness, clear image

Accurate Measurement of Dry/Wet Pulp Chamber

The new generation of core strengths improves the accuracy rate to 98%. Accurate determination of root canal length regardless of tooth shape.
Whether the root canal is dry or has electrolyte solution left in it.

Double Protection, All-Round

Enlarged display of the apical narrowing area.
Root tip stop alarm function.

R7 Accessories List

No.	Name	Quantity
1	Main unit	1
2	File clip	4
3	Measuring wire	1
4	Battery	1
5	Lip hook	4
6	File probe	4
7	ACadapter	1
8	Instruction manual	1
9	Tester	1
10	Endo ruler	1
RAL-001 Edition:V1.1 Modity20210908		



File Clip
992909
\$19.95



File Hook
992907
\$19.95



File Probe
992906
\$19.95





Cord
992908
\$19.95

Apex Locater Accessories

Type	Price	Code
File Clip	\$19.95 each	992909
File Hook	\$19.95 each	992907
File Probe	\$19.95 each	992906
Cord	\$19.95 each	992908

ENDODONTICS



ENDODONTICS

Nu²medical