

EVACUATION

EVACUATION

NuMedical



\$21.95 /bag

10pcs/bag

Aspirator Tips - Angled 16mm

Made from high quality plastic, suitable for cheek retraction. Autoclavable up to 135°. Assorted colours of blue, green, yellow, purple and lime.

Size	Code	Colour
Adult 24mm	991281	Assorted
Adult 20mm	991282	Assorted
Pedo 16mm	991283	Assorted




\$21.95 /bag

10pcs/bag

Aspirator Tips - Curved 16mm

Made from high quality plastic, suitable for cheek retraction. Autoclavable up to 135°. Assorted colours of blue, green, yellow, purple and lime.

Size	Code	Colour
Adult 22mm	991284	Assorted
Adult 20mm	991285	Assorted
Pedo 16mm	991286	Assorted




\$21.95 /bag

10pcs/bag

Aspirator Tips - Curved

Made from high quality plastic, suitable for cheek retraction. Autoclavable up to 135°. Assorted colours of blue, green, yellow, purple and lime.

Size	Code	Colour
Adult 20mm	991287	Assorted
Pedo 16mm	991288	Assorted



\$5.95 /bag

25pcs/bag

Surgical Aspirator Tips, angled 30°

Made from 100% durable polyurethane, latex-free. Fits 11mm suction hose. Autoclavable up to 135°C. Single-use.

Size	Code	Colour
1/8" tip	991248	White
1/4" tip	991249	Green
1/16" tip	991250	Blue



\$39.95 /box

1000pcs/box

Disposable HVE Aspirator Tips

Vented HVE tips designed with unique 'S' tip and non-sharp edges for superior comfort. Assorted colours of purple, yellow, lime, blue and pink. Autoclavable up to 135°C. Single-use. Size: 145mm (L) x 11mm (Diameter)

Code	Colour
991303	Assorted

FREE SAMPLES!



\$29.95 /box

1000pcs/box

Disposable HVE Aspirator Tips

Vented HVE tips designed with unique 'S' tip and non-sharp edges for superior comfort. Autoclavable up to 135°C. Size: 145mm (L) x 11mm (Diameter)

Code	Colour
991305	White

FREE SAMPLES!

EVACUATION



**\$46.95
/box**

1000pcs/box

Saliva Ejectors, Assorted

Disposable saliva ejectors for high volume evacuation of saliva and water from the oral cavity. Bendable with non-removable tip. Clear tips with assorted colour tubes of blue, green, yellow and purple. Size: 150mm (L) x 6.5mm (Diameter)

Code
991302



**\$34.90
/box**

1000pcs/box

Saliva Ejectors, Blue Tip/Clear Tube

Disposable saliva ejectors for high volume evacuation of saliva and water from the oral cavity. Bendable with non-removable tip. Size: 142mm (L) x 6.5mm (Diameter)

Code
991306



Mint flavor

**\$5.90
/bag**

1pc x100packs
/bag

Saliva Ejectors with Mint Flavor

Ideal for supplying a good smell environment to patient in dental operation. Disposable saliva ejectors for high volume evacuation of saliva and water from the oral cavity. Bendable with non-removable tip. Size: 142mm (L) x 6.5mm (Diameter)

Code
991301



**\$18.95
each**

Fraser Suction Tips

Fraser Suction Tips No.8/10/12.

Type	Code
#8	996886
#10	996887
#12	996888



**\$25.25
each**

Saliva ejectors
not included.




Saliva Ejector Dispenser,

Ideal for storing saliva ejectors or HVE tips. Size: 167mm (L) x 150mm (W) x 120mm (H)
Cleaning instructions:
Use only soap and water on acrylic surface. Alcohol or ammonia-based disinfectants may cause cracks to appear.

Code
991253

Colour
Clear



Fits HVE Tips



**\$29.95
each**

Saliva ejectors
and HVE tips
not included.

Saliva Ejector Dispenser

Ideal for storing saliva ejectors or HVE tips. Size: 170mm (H)
Cleaning instructions:
Use only soap and water on acrylic surface. Alcohol or ammonia-based disinfectants may cause cracks to appear.

Code
991244


Colour
Clear

EVACUATION

NuMedical

EVACUATION

EVACUATION



Aspirator Brushes

Stainless steel aspirator brushes with rigid bristles for effective cleaning of evacuation tubes. Autoclavable up to 135°.

Size	Length	Diameter	Price	Code
X-Fine	165mm	5mm	\$1.65 each	991256
Fine	170mm	8mm	\$1.65 each	991257
X-Fine, Long	225mm	5mm	\$1.95 each	991255
Medium	225mm	14mm	\$2.15 each	991258
Large	225mm	20mm	\$2.15 each	991259



\$9.55
/bag

72pcs/bag

Saliva Ejector Traps

Standard No: 5508. Disposable.

Code	Colour
991265	White



\$2.85
each

Snap-On Saliva Ejector Tips

Made from high quality silicone. Autoclavable up to 135°C.

Code	Colour
991276	Grey



\$21.95
each

Metal Suction Tube

Ideal for hi-volume evacuation of saliva and water from the oral cavity. Fits 11mm suction holders. Autoclavable up to 135°. Size: 130mm (L) x 11mm (Diameter)

Code	Colour
991235	Chrome



\$6.75
/bag

10pcs/bag

Suction Tip Sleeves

Made from soft silicone material Compatible with Metal Suction Tube (#991235). Autoclavable up to 135°. Size: 40mm (Diameter)



Code	Colour
991236	Blue

EVACUATION

EVACUATION

Nu
medical



\$3.90
each

12pcs/bag

Cuspidor Strainers B

Standard No: 6200. Fits Siemens Planmeca and many cuspidor with small openings. Disposable.

Code
991277



\$3.90
each

12pcs/bag

Cuspidor Strainers A

Standard No: 6100. Fits most cuspidor bowls such as A-Dec Ampco Belmont Dabi Atlante Dansereau DCI Den-Tal-Ez Forest Knight Marus MDT Proma Ritter Schein. Disposable.

Code
991268




\$3.90
/bag

12pcs/bag

Cuspidor Strainers C

Standard No: 6400. Fits: A-Dec Cascade and Performer cuspidor bowls. Disposable.

Code
991278



\$6.95
/bag

12pcs/bag

Disposable Traps E

Standard No: 5509. Fits: DCI/Belmont Reality models, System 6000 models. Size: 38mm (Diameter)

Code
991280

Colour
White



\$6.95
/bag

12pcs/bag

Disposable Traps D

Standard No: 5500. Fits: A-Dec, Cascade, Performer. Size: 31mm (Diameter)

Code
991279

Colour
White



\$6.95
/bag

12pcs/bag

Disposable Traps D

Standard No: 5500. Fits: A-Dec, Cascade, Performer. Size: 31mm (Diameter)

Code
991298

Colour
Red

EVACUATION

EVACUATION

NuMedical



Disposable Traps A

Standard No: 5501. Used to replace the existing trap in your evacuation unit or cuspidor to reduce the risks of infection and cross-contamination, while improving the performance of your evacuation system. Size: 54mm (Diameter)

Code	Colour
991263	White
991293	Blue



Disposable Traps B

Standard No: 5503. Used to replace the existing trap in your evacuation unit or cuspidor to reduce the risks of infection and cross-contamination, while improving the performance of your evacuation system. Size: 48mm (Diameter)


Code	Colour
991264	White
991294	Green



Disposable Traps C

Standard No: 5505. These disposable traps are designed to replace the existing trap in your evacuation unit or cuspidor to reduce the risks of infection and cross-contamination, while improving the performance of your evacuation system. Fits: Pelton & Crane with white plastic canisters. Size: 95mm (L)


Code	Colour
991275	White
991297	Brown



Disposable Traps F

Used to replace the existing trap in your evacuation unit or cuspidor to reduce the risks of infection and cross-contamination, while improving the performance of your evacuation system. Size: 73mm (Diameter)


Code	Colour
991289	Beige
991290	White



Disposable Traps G

Standard No: 5502. Used to replace the existing trap in your evacuation unit or cuspidor to reduce the risks of infection and cross-contamination, while improving the performance of your evacuation system. Size: 57mm (Diameter)

Code	Colour
991291	Yellow
991292	White




Disposable Traps H

Used to replace the existing trap in your evacuation unit or cuspidor to reduce the risks of infection and cross-contamination, while improving the performance of your evacuation system. Size: 70mm (Diameter)

Code	Colour
991295	Grey
991296	White

EVACUATION

EVACUATION




\$15.65
/bag

12pcs/bag

Cuspidor Strainers, 1

Fits Adec, Cascade and Performer.
Disposable.

Code	Colour
991300	White



\$15.65
/bag

12pcs/bag

Cuspidor Strainers, 2

Fits: A-Dec, Ampco, Belmont, Dab
Atante, Danserau, Forest, Knight,
Marus, MDT, Proma, Ritter and Schein.
Disposable.

Code	Colour
991299	White



\$81.90
/box

12pcs/box

Disposable Canisters A

Filters for central vacuum pumps, threads-internal. Fits: Air Technigues-Vacstar Models, Apollo 1 hp, 1.5 Single: Tandem Models, 1, 2hp Den-Tal-EZ Single Models, Ohmeda-Matrix most models, Turbine Industries, Techwest-Whirl Vac Models, Vacudent-Workhorse Single and other units.

Code	Colour
991266	Clear



\$81.90
/box

8pcs/box

Disposable Canisters B

Filters for central vacuum pumps, threads-Internal, Fits: Den-Tal-Ez dual models made after 1978, Dentsply made after 1991, Matrix 1000, 2000 and Drymax, Vortex Vacuum Pump (JDS).

Code	Colour
991267	Clear



\$9.95
each

HVE Lever, 11mm

Disposable HVE lever. Cleaning instruction: Clean with disinfectant or alcohol wipes.

Code
991239



\$9.95
each

HVE Lever, 16mm

Disposable HVE lever. Cleaning instruction: Clean with disinfectant or alcohol wipes.

Code
991238

Nu
medical

EVACUATION

EVACUATION

Nu2medical



\$19.00
each

Rotary Saliva Ejector Valve

Autoclavable up to 135°.

Code	Colour
991272	Grey




\$19.00
each

Rotary HVE Lever

Autoclavable up to 135°.

Code	Colour
991270	Grey



\$19.00
each

HVE Valve Lever

Autoclavable up to 135°.

Code	Colour
991269	Grey



\$19.00
each

Saliva Ejector Valve Lever

Autoclavable up to 135°.

Code	Colour
991271	Grey



\$2.95
each

HVE Lever

Autoclavable up to 135°.

Code	Colour
991273	Grey



\$2.95
each

Saliva Ejector Lever

Autoclavable up to 135°.

Code	Colour
991274	Grey



\$5.95
/bag

5pcs/bag

Suction Tube Adaptors A

Adaptor to connect 11mm suction tubes to European HVE 16mm hose. Autoclavable up to 135°C.

Code	Colour
991260	White



\$5.00
each

Suction Tube Adaptors

Adaptor to connect saliva ejectors (6.5mm) to standard HVE hose 11mm. Autoclavable up to 135°C.

Code	Colour
991261	White

EVACUATION



\$9.75
/bag

50pcs/bag

Saliva Ejector Cushions

These cushions are used together with saliva ejectors. The cushion fits and slides easily over the saliva ejector to help prevent tissue grabbing or irritation lowers suction noise and provides efficient suction and evacuation with maximum patient comfort. Made from non-latex open cell foam. Size: 30mm (L) x 15.8mm (W).

Code
990905



\$26.65
/box

200pcs/box

Mirror Wipes

Remove moisture and debris from the mouth mirror easily and quickly. With these handy mirror wipes dental professionals do not need to remove the mirror from the patient's mouth. These mirror wipes can be affixed easily onto syringe tips.

Code
990909



\$59.95
each

Outer diameter: 21mm
Inner diameter: 17mm



Nu Hve Suction Tube Grey

Ideal for HVE Suction. Soft feeling. Internal Dia.17 mm and External Dia.21 mm. The length arrives 1.7meters.

Code
991307

Colour
Grey




Saliva Ejector Holoder

Ideal for ultrasonic scaling and sealants! Simply slide holder onto your saliva ejector to provide continuous, comfortable, hands-free suction with no additional attachments required. By placing buccal side, the holder confidently holds the saliva ejector right where you need it – the back corner of the patients mouth. Easy placement, Patient comfort, Single-use, Latex-Free.

\$13.95
/box

50pcs/box



Code
992763

EVACUATION

NuMedical