

GLOVES AND DENTAL DAMS

GLOVES AND DENTAL DAMS




Scan it for price
300pcs/box

Nusuperslim Cobalt Blue Powder Free Nitrile Examination Glove

Soft like skin. Cobalt blue powder free nitrile examination glove. Non-sterile. Ambidextrous, finger textured and beaded cuff. For single use only. 300pcs/box. The length is almost 240mm which is longer than nitrile gloves in the market.

FREE SAMPLES!

Size	Code	Colour
XS	990040	Blue
S	990041	
M	990042	
L	990043	



Scan it for price
100pcs/box

NuEco Latex Exam Gloves, Powder Free

Economical online-chlorinated powder free latex exam gloves by NuMedical!

FREE SAMPLES!

Size	Code	Colour
XS	990019	Natural Latex
S	990020	
M	990021	
L	990022	



Scan it for price
200pcs/box

FREE SAMPLES!

Basic Nitrile Examination Glove

Soft and stretchy nitrile formulation with excellent tactile sensitivity. Latex-free

Size	Code	Colour
XS	990060	Blue
S	990061	
M	990062	
L	990063	



Scan it for price
100pcs/box


FREE SAMPLES!

Black Lion Nitrile Examination Glove

Soft and stretchy nitrile formulation with excellent tactile sensitivity. Latex-free

Size	Code	Colour
S	990027	Black
M	990028	
L	990029	

NuMedical



\$29.65 each



Fits box dimensions:
≤ 235mm (L)


Gloves and Tissue Box Holder Type 1

Ideal for storing gloves or tissue. Pre-drilled holes for wall-mounting. Size: 250mm (L) x 125mm (W) x 110mm (H)

Cleaning instructions:

Use only soap and water on acrylic surface. Alcohol or ammonia-based disinfectants may cause cracks to appear.

Code	Colour
990005	Clear



\$29.65 each



Fits box dimensions:
≤ 235mm (L)

Gloves and Tissue Box Holder Type 2

Ideal for storing gloves or tissue. Pre-drilled holes for wall-mounting. Size: 250mm (L) x 130mm (W) x 90mm (H)

Cleaning instructions:


Use only soap and water on acrylic surface. Alcohol or ammonia-based disinfectants may cause cracks to appear.

Code	Colour
990006	Clear

GLOVES AND DENTAL DAMS

GLOVES AND DENTAL DAMS

NuMedical



\$37.35
each



**Fits box
dimensions:
≤ 235mm (L)**


Gloves and Tissue Box Holder Type 3

Ideal for storing gloves or tissue. Pre-drilled holes for wall-mounting. Size: 250mm (L) x 90mm (W) x 250mm (H)

Cleaning instructions:

Use only soap and water on acrylic surface. Alcohol or ammonia-based disinfectants may cause cracks to appear.

Code	Colour
990007	Clear



\$47.25
each



**Fits box
dimensions:
≤ 235mm (L)**


Gloves and Tissue Box Holder Type 4

Ideal for storing gloves or tissue. Pre-drilled holes for wall-mounting. Size: 250mm (L) x 90mm (W) x 370mm (H)

Cleaning instructions:

Use only soap and water on acrylic surface. Alcohol or ammonia-based disinfectants may cause cracks to appear.

Code	Colour
990008	Clear



\$18.55
each



**Fits box
dimensions:
≤ 135mm (H)**


Gloves Box Holder Type 5

Ideal for storing a single gloves box. Pre-drilled holes for wall-mounting. Size: 250mm (L) x 102mm (W) x 145mm (H)

Cleaning instructions:

Use only soap and water on acrylic surface. Alcohol or ammonia-based disinfectants may cause cracks to appear.

Code	Colour
990009	Clear



\$23.95
each



**Fits box
dimensions:
≤ 132mm (H)**

Gloves Box Holder Type 6

Ideal for storing two gloves boxes. Pre-drilled holes for wall-mounting. Size: 250mm (L) x 102mm (W) x 275mm (H)

Cleaning instructions:

Use only soap and water on acrylic surface. Alcohol or ammonia-based disinfectants may cause cracks to appear.

Code	Colour
990010	Clear



\$25.65
each



**Fits box
dimensions:
≤ 131mm (H)**

Gloves Box Holder Type 7

Ideal for storing three gloves boxes. Pre-drilled holes for wall-mounting. Size: 250mm (L) x 102mm (W) x 405mm (H)

Cleaning instructions:

Use only soap and water on acrylic surface. Alcohol or ammonia-based disinfectants may cause cracks to appear.

Code	Colour
990011	Clear

GLOVES AND DENTAL DAMS

GLOVES AND DENTAL DAMS

NuMedical

\$4.35 each

Rubber Dam Molar Clamps (WITHOUT Dam-Engaging Projections)

					
#15 996500	#18 996501	#24 996502	#25 996503	#26 996504	#28 996505

				
#30 996506	#31 996507	#51 996508	#138 996509	#139 996510

Rubber Dam Molar Clamps, B4 For Anterior Teeth


#B4 996878

Rubber Dam Molar Clamps (WITHOUT Dam-Engaging Projections)

				
#27N 996856	#8A 996857	#14A 996858	#12A 996859	#13A 996860

Rubber Dam Molar Clamps \$7.95 each


#44(#B4 GGold) 997851

Rubber Dam Molar Clamps (WITH Dam-Engaging Projections)

					
#200 996511	#201 996512	#202 996513	#203 996514	#204 996515	#205 996543

Rubber Dam Biscupid Clamps (WITHOUT Dam-Engaging Projections)

		
#22 996523	#27 996524	#29 996525

Rubber Dam Labial Clamps (WITH Dam-Engaging Projections)

	
#210 996520	#211 996521

Rubber Dam Biscupid Clamps (WITH Dam-Engaging Projections)

			
#206 996516	#207 996517	#208 996518	#209 996519

Rubber Dam Cervical Clamp


#212 996522

GLOVES AND DENTAL DAMS

GLOVES AND DENTAL DAMS

Model	Quantity	Photo	Tooth Type	Function	Tooth Position
H1-1	1		Molar with arms A	Suitable for upper and lower molars, and arms could fix rubber dam during multi-tooth treatment	
H1-2	1		Molar with arms B	Suitable for upper and lower molars, and arms could fix rubber dam during multi-tooth treatment	
H1	3		Molar	Suitable for upper and lower molars, especially for crowns with complete shape	
T1	3		Bicuspid	Suitable for upper and lower bicuspids, no wings, with fixing band to fix rubber dam during treatment	



\$79.95
/box

UDG Concept One 3D Rubber Dam Clamps Kit

Ideal for upper and lower molars, and arms could fix rubber dam during multi-tooth treatment.

Code
997861



\$16.85
/box

52pcs/box

Cranberry Dental Dams, 5" x 5"

Eliminate the interference of saliva to treatment.
Clear vision to improve work efficiency.
Strong moisture barrier to prevent infection and accidental inhalation.
Size: side length 127mm.
Class I medical devices; Single-use.
High temperature and high pressure sterilization is not allowed.

Code	Scent
991051	Plain Scent
991053	Mint



\$16.85
/box

36pcs/box

Cranberry Dental Dams, 6" x 6"

Eliminate the interference of saliva to treatment.
Clear vision to improve work efficiency.
Strong moisture barrier to prevent infection and accidental inhalation.
Size: side length 152mm.
Class I medical devices; Single-use.
High temperature and high pressure sterilization is not allowed.

Code	Scent
991052	Plain Scent
991054	Mint



\$25.95
/box

15pcs/box

Non-Latex Spearmint Dental Dam, 6" x 6"

Perfect cut for an uniformed thickness and strength to provide consistent workability. Unparalleled tensile strength for uncompromising tear resistance >600%. High contrast for optimal operation visual. Excellent recovery properties for effective isolation.

Code
990970



\$6.75
each

Dental Dam Frame - Plastic

Unique design with strong grip and enhanced patient comfort.
Autoclavable up to 135 degree. Reusable.

Code
990998

Nu
medical

GLOVES AND DENTAL DAMS


GLOVES AND DENTAL DAMS



Gentle for sensitive dental neck,
ceramic teeth and all porcelain crowns



Anti-slip processed in the area
contacting with tooth surface



Premolar

Molar

Molar

\$39.95
/box


6pcs/box

Resin Dental Dam Clamps

2 Premolars, 2 Molars, 2 Molars

To protect sensitive tooth, versatile design fits a
variety of crowns and easy to retain the dental dam
With good elasticity, resilience and retention
Gentle for sensitive dental neck, ceramic teeth and all
porcelain crowns

Code
998016



Place the clamp on the premolar with dental
clamp forceps



Short

Long

\$25.95
/box

4pcs/box

Distal-Extension Dental Dam Clamp

2 long, 2 Shorts

For fixing the dental dam in cases of distal-extension
class II restoration

Place the clamp on the premolar with dental dam
clamp on the premolar with dental dam clamp forceps
When using the dental dam on the distal-extension
the clamp fixes the dental dam and helps to leave
enough treatment area

Code
998017



Fully expose the distal treatment area,
not conflicting with the instruments



\$7.95
/box

2.1m/box

Dental Dam Stabilizing Cord

Non-latex, with good resilience, to
stabilize the dental dam

Code
998019

Nu2medical

GLOVES AND DENTAL DAMS

GLOVES AND DENTAL DAMS

Numedical




\$9.95
each

Dental Dam Frame - Stainless Steel

Made from high grade stainless steel.
Autoclavable up to 135 degree. Reusable.
Size: 100mm(L) x 100mm(W).

Code
990997



\$9.95
each

Dental Dam Frame - Stainless Steel

Made from high grade stainless steel.
Autoclavable up to 135 degree. Reusable.
Size: 127mm (L) x 127mm (W)

Code
991000



\$42.95
each

Rubber Dam Punch, Ainsworth

Ainsworth-style rubber dam punch.
Size: 160mm (L)

Code
996762



\$42.95
each

Rubber Dam Punch, Ivory

Ivory-style rubber dam punch.
Size: 170mm (L)

Code
996763




\$35.95
each

Rubber Dam Forceps

Brewer-style rubber dam punch.
Size: 175mm (L)

Code
997852



\$35.95
each

Rubber Dam Forceps, Ivory Lightweight

Ivory-style rubber dam punch.
Size: 170mm (L)

Code
996765