

INSTRUMENTS

PERSONALISE Your Own Stainless Steel Mouth Mirrors

1 SELECT MIRROR TYPE



2 SELECT HANDLE



1 SELECT MIRROR TYPE

	Size	Quantity	Price	Code
Ultra Vision Front Surface, No. 4	22mm	12pcs/bag	\$29.95/bag	996778
Ultra Vision Front Surface, No. 5	24mm	12pcs/bag	\$29.95/bag	996779
Plain Face, No. 3	20mm	12pcs/bag	\$14.95/bag	996784
Plain Face, No. 4	22mm	12pcs/bag	\$14.95/bag	996785
Plain Face, No. 5	24mm	12pcs/bag	\$14.95/bag	996786
Magnifying, No. 4	22mm	12pcs/bag	\$19.95/bag	996758
Magnifying, No. 5	24mm	12pcs/bag	\$19.95/bag	996759
Front Surface, No. 4	22mm	12pcs/bag	\$19.95/bag	996760
Front Surface, No. 5	24mm	12pcs/bag	\$19.95/bag	996761

2 SELECT HANDLE

	Price	Code
Mirror Handle, Round. 120mm (L) x 8mm (Diameter)	\$3.95 each	996780
Mirror Handle, Octagonal. 120mm (L) x 8mm (Diameter)	\$3.95 each	996781
Mirror Handle with Ruler. 120mm (L) x 8mm (Diameter)	\$3.95 each	996782
Mirror Handle, Extra Light. 120mm (L) x 8mm (Diameter)	\$6.95 each	996783
Mirror Handle, XL. 120mm (L) x 8mm (Diameter)	\$12.95 each	996787
Mirror Handle, XL with titanium-coated stem. 120mm (L) x 8mm (Diameter)	\$17.95 each	996788

INSTRUMENTS

PERSONALISE Your Own Stainless Steel Mouth Mirrors

1

SELECT MIRROR TYPE



2

SELECT HANDLE



1

SELECT MIRROR TYPE

	Size	Quantity	Price	Code
Ultra Vision Front Surface, No. 4	22mm	12pcs/bag	\$29.95/bag	996823
Ultra Vision Front Surface, No. 5	24mm	12pcs/bag	\$29.95/bag	996824
Plain Face, No. 3	20mm	12pcs/bag	\$14.95/bag	996825
Plain Face, No. 4	22mm	12pcs/bag	\$14.95/bag	996826
Plain Face, No. 5	24mm	12pcs/bag	\$14.95/bag	996827
Magnifying, No. 4	22mm	12pcs/bag	\$19.95/bag	996828
Magnifying, No. 5	24mm	12pcs/bag	\$19.95/bag	996829
Front Surface, No. 4	22mm	12pcs/bag	\$19.95/bag	996830
Front Surface, No. 5	24mm	12pcs/bag	\$19.95/bag	996831

2

SELECT HANDLE

	Price	Code
Mirror Handle, Round. 120mm (L) x 8mm (Diameter) 	\$3.95 each	996832
Mirror Handle, Octagonal. 120mm (L) x 8mm (Diameter) 	\$3.95 each	996833
Mirror Handle with Ruler. 120mm (L) x 8mm (Diameter) 	\$3.95 each	996834
Mirror Handle, Extra Light. 120mm (L) x 8mm (Diameter) 	\$6.55 each	996835
Mirror Handle, XL. 120mm (L) x 8mm (Diameter) 	\$12.95 each	996836
Mirror Handle, XL with titanium-coated stem. 120mm (L) x 8mm (Diameter) 	\$17.95 each	996837

INSTRUMENTS

NuMedical

INSTRUMENTS

INSTRUMENTS

NuMedical



\$79.95
/box

10pcs/box



Mouth Mirror Anit-Fog Simple Stem, Flat Handle

Mouth Mirror Anti-Fog Simple Stem Type, Ideal for supplying good visual in dental procedure. Compatible with stem type handle, like 996832 - Round, 996833 - Octagonal, 996834 - With Rule, 996835 - Extra Light, 996836 - XL, 996837 - XL with Titanium-coated.

Code
997857



\$79.95
/box

10pcs/box



Mouth Mirror Anit-Fog Simple Stem, Round Handle

Mouth Mirror Anti-Fog Simple Stem Type, Ideal for supplying good visual in dental procedure. Compatible with stem type handle, like 996832 - Round, 996833 - Octagonal, 996834 - With Rule, 996835 - Extra Light, 996836 - XL, 996837 - XL with Titanium-coated.

Code
997858



\$19.95
each



Plaster Forceps

Used for handling and cutting plaster models.

Code
996746



\$19.95
each

Amalgam Well

Stainless steel amalgam well with anti-skid silicone "O" ring base.

Code
996741



\$14.95
each

Measuring Gauge, Boley

Used to measure tooth dimension (length, width, thickness) in millimeter increments.

Code
996739



\$14.95
each

Measuring Gauge, Iwanson

Used to measure tooth dimension (length, width, thickness) in millimeter increments.

Code
996740

INSTRUMENTS



\$9.95
each

Scissors, Stainless Steel, Curved

Scissors, Stainless Steel, Curved Festooning for 3m Espe 801203 Users.

Code
996877



\$68.95
each

Needle Holder Castroviejo with Tungsten Carbide Straight 14.5cm

Needle Holder Castroviejo with Tungsten Carbide, Straight 14.5cm


Code
996899



\$14.95
each

Bolten Surgical Towel Clips

Code
996884



\$29.95
each

Crown and Bridge Remover

For removal of crowns, bridges and inlays. Comes with 3 interchangeable tips.

Code
996752



\$2.95
each

Scalpel Handle #3 with Ruler

Stainless Steel scalpel handle #3 with ruler, 0 - 6cm measurement.

Code
996558



\$29.95
each

Bone Rongeur, Friedman, 140mm

Stainless Steel bone rongeur, Friedman. Length: 140mm

Code
996559

INSTRUMENTS

NuMedical

INSTRUMENTS

INSTRUMENTS



\$4.95
each



Scissors, Iris Straight

Stainless Steel iris straight scissors.
Length: 110mm

Code
996544



\$4.95
each



Scissors, Iris Curved

Stainless Steel iris curved scissors.
Length: 110mm

Code
996545



\$29.95
each



Scissors, Iris Straight

Tungsten Carbide iris straight scissors.
Length: 115mm

Code
996546



\$29.95
each



Scissors, Iris Curved

Tungsten Carbide iris curved scissors.
Length: 115mm

Code
996547



\$39.95
each



Suture Needle Holder, Hegar

Tungsten Carbide suture needle holder, Hegar.
Length: 140mm

Code
996548



\$39.95
each



Suture Needle Holder, Mathieu

Tungsten Carbide suture needle holder, Mathieu.
Length: 140mm

Code
996549

NuMedical


INSTRUMENTS

INSTRUMENTS

NuMedical



\$13.95
each



Suture Needle Holder, 135mm

Stainless Steel suture needle holder.
Length: 135mm

Code
996550



\$13.95
each



Suture Needle Holder, 150mm

Stainless Steel suture needle holder.
Length: 135mm

Code
996551



\$13.95
each



Splinter Forceps, Straight, 140mm

Stainless Steel splinter forceps, straight.
Length: 140mm

Code
996560



\$13.95
each



Splinter Forceps, Curved, 140mm

Stainless Steel splinter forceps, curved.
Length: 140mm

Code
996561



\$9.95
each



Hemostatic Forceps, Straight, 140mm

Stainless Steel hemostatic forceps, straight.
Length: 140mm

Code
996555



\$9.95
each



Hemostatic Forceps, Curved, 140mm

Stainless Steel hemostatic forceps, curved.
Length: 140mm

Code
996556

INSTRUMENTS

INSTRUMENTS

NuMedical



\$9.95
each



Hemostatic Forceps, Straight, 120mm

Stainless Steel hemostatic forceps, straight.
Length: 120mm

Code
996553



\$9.95
each



Hemostatic Forceps, Curved, 120mm

Stainless Steel hemostatic forceps, curved.
Length: 120mm

Code
996554



\$13.95
each



Suture Needle Holder, 140mm

Stainless Steel suture needle holder.
Length: 135mm

Code
996552




\$12.95
each

Periosteal Elevator Molt, No.9

The rounded edge of the periosteal instrument minimises tissue tearing when lifting the gum flap while the sharper end makes the initial lifting of the flap easier.

Code
996557



\$5.95
each

Adson Soft Tissue Forceps

Code
996876



\$14.95
each

Atraumatic Tweezers/Forceps Flat Handle

Code
996893

INSTRUMENTS

INSTRUMENTS

NuMedical



\$19.95
each

Plaster Knife, 5R

Stainless steel knife with wooden handle used for handling, shaping and cutting plaster model. Size: 160mm (L)

Code
996747



\$19.95
each

Plaster Knife, 5R

Stainless steel knife with wooden handle used for handling, shaping and cutting plaster model. Size: 160mm (L)

Code
996749



\$9.95
each

Gillies Tweezer

Tissue Tweezers, Gillies 15cm

Code
996882



\$5.95
each

Wax Knife, Fahren Small

Stainless steel knife with wooden handle used for wax carving and mixing.

Code
996774



\$19.95
each

Plaster Knife, 5R

Stainless steel knife with wooden handle used for handling, shaping and cutting plaster model. Size: 160mm (L)

Code
996748



\$5.95
each

Wax Knife, Fahren Large

Stainless steel knife with wooden handle used for wax carving and mixing.

Code
996775

INSTRUMENTS

INSTRUMENTS



\$73.45
each

Amalgam Capsule Plier

Ideal for fitting GC & 3M Capsule.

Code
993149



\$28.55
each

Tongue Retractor

Designed to hold patient's tongue away from the surgical area. Autoclavable up to 135°C.
Size: 190mm (L)

Code
993265

3mm



4mm



5mm



6mm



Sterile Disposable Biopsy Punch

NuMedical's sterile disposable biopsy punches have a seamless, stainless steel and razor sharp cutting edge. They are individually packaged in a sterilized package. Its ribbed handle provides maximum control and comfort.

\$49.50
/box

10pcs/box

Code	Size	Color
996889	3mm	White
996890	4mm	Blue
996891	5mm	Pink
996892	6mm	Green

\$3.95
each



Wax Carver, Beale Large

Used in wax carving and mixing.

Code
996773

\$3.95
each



Wax Carver, Beale Small

Used in wax carving and mixing.

Code
996772

NuMedical

INSTRUMENTS

INSTRUMENTS



\$3.95
each

Wax Carver, Le Cron

Used in wax carving and mixing.

Code
996769



\$3.95
each

Wax Carver, Zahle

Used in wax carving and mixing.

Code
996770




\$3.95
each

Wax Carver, 5

Used in wax carving and mixing.


Code
996771

Nu2medical



INSTRUMENTS

INSTRUMENTS



NuMedical

\$21.95
each

Extraction Forceps, Pedo

<p>ANATOMICAL HANDLE</p> <p>Anatomical Handle</p>	<p>Upper Anterior</p> <p>996598</p>	<p>Upper Roots Type 1</p> <p>996599</p>	<p>Upper Molars</p> <p>996600</p>	<p>Upper Roots Type 2</p> <p>996601</p>
	<p>Lower Anterior</p> <p>996602</p>	<p>Lower Molars</p> <p>996603</p>	<p>Lower Roots</p> <p>996604</p>	

INSTRUMENTS

Extraction Forceps

 <p>Upper Incisors and Bicuspid, Type 1</p> <p>996605 \$21.95</p>	 <p>Upper Incisors and Bicuspid, Type 2</p> <p>996606 \$21.95</p>	 <p>Upper Incisors and Bicuspid, Type 3</p> <p>996607 \$21.95</p>	 <p>Upper Incisors and Bicuspid, Type 4</p> <p>996608 \$21.95</p>
 <p>Upper Premolars</p> <p>996609 \$21.95</p>	 <p>Lower Premolars</p> <p>996610 \$21.95</p>	 <p>Upper Molars Right</p> <p>996611 \$21.95</p>	 <p>Upper Molars Left</p> <p>996612 \$21.95</p>
 <p>Upper Molars Type 1</p> <p>996613 \$21.95</p>	 <p>Lower Molars Type 1</p> <p>996614 \$29.95</p>	 <p>Lower Molars Right</p> <p>996615 \$29.95</p>	 <p>Lower Molars Left</p> <p>996616 \$29.95</p>
 <p>Upper Incisors, Bicuspid and Roots</p> <p>996617 \$21.95</p>	 <p>Upper Roots Type 1</p> <p>996618 \$21.95</p>	 <p>Lower Roots Type 1</p> <p>996619 \$21.95</p>	 <p>Lower Roots Type 2</p> <p>996620 \$21.95</p>

INSTRUMENTS

INSTRUMENTS

NuMedical

 <p>Upper Roots Type 2</p> <p>996621 \$21.95</p>	 <p>Lower Roots Type 1</p> <p>996622 \$21.95</p>	 <p>Upper Roots Type 3</p> <p>996623 \$21.95</p>	 <p>Upper Roots Type 4</p> <p>996624 \$21.95</p>
 <p>Upper Roots Type 5</p> <p>996625 \$21.95</p>	 <p>Lower Roots Type 1</p> <p>996626 \$21.95</p>	 <p>Lower Molars Type 2</p> <p>996627 \$21.95</p>	 <p>Lower Roots incisors & Bicuspid</p> <p>996628 \$21.95</p>
 <p>Lower Roots incisors & Bicuspid</p> <p>996629 \$21.95</p>	 <p>Upper Roots Type 6</p> <p>996630 \$21.95</p>	 <p>Lower Wisdom Type 1</p> <p>996631 \$21.95</p>	 <p>Lower Wisdom Type 2</p> <p>996632 \$21.95</p>
 <p>Lower Anterior Type 3</p> <p>996633 \$29.95</p>	 <p>Lower Anterior Type 4</p> <p>996634 \$29.95</p>	 <p>Lower Anterior Type 5</p> <p>996635 \$29.95</p>	 <p>Lower Molars Type 6</p> <p>996636 \$29.95</p>

INSTRUMENTS

INSTRUMENTS


NuMedical



Lower Molars

996868

\$21.95



Lower Molars
Type 6

996636

\$29.95



Upper Roots and
Premolars

996637

\$21.95



Lower Molars

996867

\$29.95



Upper Molars Right

996869

\$21.95



Upper Molars Left

996870

\$21.95



Lengthened Upper Pre-
molar

996866

\$29.95



Lengthened Upper Roots

996865

\$29.95



Upper Premolar

996863


\$29.95



Upper Incisors

996862

\$29.95



Lower Molars
#22

996864


\$21.95

INSTRUMENTS

INSTRUMENTS

Root Elevator, Bein

\$11.95
each



Type 1	Type 2	Type 3	Type 4	Type 5	Type 6
996562	996563	996564	996565	996566	996567

Root Elevator, Cryer

\$11.95
each

Elevator, Straight

\$11.95
each



Type 1	Type 2	Type 3	Type 4
996568	996569	996570	996571



#301
996871

Root Elevator, Apical

\$11.95
each



Type 1	Type 2	Type 3	Type 4
996578	996579	996580	996581

Nu2medical

INSTRUMENTS

INSTRUMENTS

Nu²medical

Root Elevator, Warwick James

\$11.95 each



Root Elevator, Coupland

\$11.95 each





Root Elevator, Heidbrink B

\$11.95 each



Root Elevator, Flohr

\$21.95 each



Root Elevator, Luxating Type

\$21.95 each




INSTRUMENTS

INSTRUMENTS

Root Elevator, Heidbrink A

\$11.95
each



NuMedical



Right

Elevator, Potts

Code	Size
996872	Right

\$21.95
each



Left

Elevator, Potts

Code	Size
996873	Left


















\$21.95
each

INSTRUMENTS

Excavators, Double-Ended, Hollow Handle

\$4.95 each

Code

Type 1		996643
Type 2		996644
Type 3		996645
Type 4		996646
Type 5		996647
Type 6		996648
Type 7		996649
Type 8		996650
Type 9		996651
Type 10		996652
Type 11		996653
Type 12		996654
Type 13		996655
Type 14		996656
Type 15		996657
Type 16		996658
Type 17		996659

INSTRUMENTS

NuMedical

INSTRUMENTS






INSTRUMENTS

NuMedical

Probes, Hollow Handle

\$10.20 each






Code

Probe, Hollow Handle, Type 1		996638
Probe, Hollow Handle, Type 2		996639
Probe, Hollow Handle, Type 3		996640
Probe, Hollow Handle, Type 4		996641
Probe, Hollow Handle, Type 5		996642

Probes

\$5.45 each

Code

Probe, Type 1		996661
Probe, Type 2		996662
Probe, Type 3		996663
Probe, Type 4		996664
Probe, Type 5		996665

Probe, Double-Ended

\$5.95 each

Code

Probe, Double-Ended		996660
---------------------	--	--------

Explorer, Double-Ended

\$10.95 each

Code

Explorer, Double-Ended		996542
------------------------	--	--------

Endo Plugger

\$4.95 each

Code

Endo Plugger		996844
--------------	--	--------

INSTRUMENTS

INSTRUMENTS

Perio Probes with Laser Marking, Double-Ended

\$11.95 each

Code

Perio Probe, Double-Ended,
23/Williams, 10/9/8/7/5/3/2/1



996682

Perio Probe, Double-Ended,
23/WHO, 11.5/8.5/5.5/3.5



996683

Perio Probe, Double-Ended,
23/11, 11/8/6/3



996684

Perio Probe, Double-Ended,
23/12, 12/9/6/3



996685

Perio Probes with Laser Marking

\$11.95 each

Code

Perio Probe, Williams,
1/2/3/5/7/8/9/10



996686

Perio Probe, WHO,
3.5/5.5/8.5/11.5



996687

Perio Probe,
11, 3/6/8/11



996688

Perio Probe,
12, 3/6/9/12



996689

Tweezers

\$7.95 each

Code

Perry Tweezer, 125mm



996789

College Tweezer, 150mm



996790

Self-Lock Tweezer, 150mm



996791

Dressing Tweezer, 145mm



996792
















NuMedical

INSTRUMENTS

Curettes and Scalers, Hollow Handle

\$11.95 each

Code

Curette, Hollow Handle, Gracey, 1-2		996526
Curette, Hollow Handle, Gracey, 3-4		996527
Curette, Hollow Handle, Gracey, 5-6		996528
Curette, Hollow Handle, Gracey, 7-8		996529
Curette, Hollow Handle, Gracey, 9-10		996530
Curette, Hollow Handle, Gracey, 11-12		996531
Curette, Hollow Handle, Gracey, 13-14		996532
Curette, Hollow Handle, Columbia, 2L-2R		996533
Curette, Hollow Handle, Columbia, 4L-4R		996534
Curette, Hollow Handle, McCall, 13S-14S		996535
Scaler, Hollow Handle, Sickle, Type 1		996536
Scaler, Hollow Handle, Sickle, Type 2		996537
Scaler, Hollow Handle, Sickle, Type 3		996538
Scaler, Hollow Handle, Sickle, Type 4		996539
Scaler, Hollow Handle, Sickle, Type 5		996540

INSTRUMENTS

Spoon Excavators, Double-Ended

\$5.95 each

Code

Excavator, Double-Ended,
17, 1.2mm - 1.5mm



996677

Excavator, Double-Ended,
18, 1.5mm - 1.8mm



996678

Excavator, Double-Ended,
19, 2mm - 2.2mm



996679

Excavator, Double-Ended,
20, 2.2mm - 2.5mm



996680

Excavators, Double-Ended, Hollow Handle

\$11.95 each

Code

Excavator, Double-Ended,
Hollow Handle, 17



996708

Excavator, Double-Ended,
Hollow Handle, 18



996709

Excavator, Double-Ended,
Hollow Handle, 19



996710

Excavator, Double-Ended,
Hollow Handle, 127/128



996711

Excavator, Double-Ended,
Hollow Handle, 129/130



996712

Excavator, Double-Ended,
Hollow Handle, 153/154



996713

Excavator, Double-Ended,
Hollow Handle, 218/219



996714

Excavator, Double-Ended,
Hollow Handle, 220/221



996715

INSTRUMENTS

NuMedical

INSTRUMENTS












INSTRUMENTS

NuMedical

Amalgam Instruments, Double-Ended, Hollow Handle

\$11.95 each

Code

Plugger, Hollow Handle, Serrated, 0&1, 1.5mm - 2mm		1.5mm 996690
Plugger, Hollow Handle, Serrated, 1&2, 2mm - 2.5mm		2.0mm 996691
Plugger, Hollow Handle, Smooth, 0&1, 1.5mm - 2mm		2.0mm 996692
Plugger, Hollow Handle, Smooth, 1&2, 2mm - 2.5mm		2.0mm 996693
Burnisher Ball, 18		1.5mm 996694
Carver, Hollenback, 3S		1.0mm 996695
Carver, Hollenback, 1/2		0.75mm 996696
Carver, Interproximal, 3S		0.8mm 996697
Burnisher, Michigan, 1		1.2mm 996698
Burnisher, Michigan, 2		2.0mm 996699
Plugger, Smooth, 1&2, 1.5mm - 2.5mm		3.0mm 996841

Endo Probe

\$11.95 each

Code

General endodontic explorer for endodont applications.		996843
--	--	--------

Periosteal Elevator Buser, No.6,

\$12.95 each

Code

Periosteal Elevator Buser, No.6		4.6mm 996894
---------------------------------	--	-----------------

INSTRUMENTS

Amalgam Instruments, Double-Ended

\$5.95 each

Code

Plugger, Double-Ended,
Serrated,
0&1, 1.5mm - 2mm



2.0mm
996666

Plugger, Double-Ended,
Serrated,
1&2, 2mm - 2.5mm



2.5mm
996667

Plugger, Double-Ended,
Smooth,
0&1, 1.5mm - 2mm



2.0mm
996668

Plugger, Double-Ended,
Smooth,
1&2, 2mm - 2.5mm



2.0mm
996669

Plugger, Hollenback,
Double-Ended, 2,
2mm-3mm



2.0mm
996670

Burnisher Ball,
Double-Ended,
28&29



7.0mm
996671

Burnisher Ball,
Double-Ended,
1, 1.2mm - 1.5mm



1.2mm
996672

Burnisher Ball,
Double-Ended,
1, 2mm - 2.5mm



2.0mm
996673

Carver, Hollenback,
Double-Ended, 3




1.2mm
996674
2.4mm

Carver, Hollenback,
Double-Ended, 1/2



0.85mm
996675
1.7mm

Ward, Double-Ended, 1



1.1mm
996676
2.2mm

Amalgam Carrier

\$19.95 each

Code

Large/Regular,
3mm/2mm



996742

Cord Packer 113

\$11.95 each

Code

Cord Packer 113



996897

INSTRUMENTS

Nu
medical

INSTRUMENTS








INSTRUMENTS

NuMedical

Filling Instruments, Double-Ended, Hollow Handle

\$11.95 each









Code

Filling Instrument, Double-Ended, Hollow Handle, Type 1		2.2mm 2.2mm	996700
Filling Instrument, Double-Ended, Hollow Handle, Type 2		2.5mm 2.5mm	996701
Filling Instrument, Double-Ended, Hollow Handle, Type 3		2.0mm 2.0mm	996702
Filling Instrument, Double-Ended, Hollow Handle, Type 4		2.0mm 2.2mm 1.1mm	996703
Filling Instrument, Double-Ended, Hollow Handle, Type 5		2.4mm 2.0mm 1.0mm	996704
Filling Instrument, Double-Ended, Hollow Handle, Type 6		0.7mm 1.3mm 2.0mm 1.2mm	996705
Filling Instrument, Double-Ended, Hollow Handle, Type 7		0.7mm 1.0mm 0.7mm 1.4mm	996706
Filling Instrument, Double-Ended, Hollow Handle, Type 8		1.25mm 2.5mm 3.0mm 1.5mm	996707

Filling Instruments, D/Ended, Hollow Handle, Teflon Coating

\$23.95 each

Code

Filling Instrument, Hollow Handle, Blue Teflon Coating, Type 1		3.0mm 3.0mm	996795
Filling Instrument, Hollow Handle, Blue Teflon Coating, Type 2		2.6mm 2.6mm	996796
Filling Instrument, Hollow Handle, Blue Teflon Coating, Type 3		2.0mm 1.9mm	996797
Filling Instrument, Hollow Handle, Blue Teflon Coating, Type 4		2.2mm 1.7mm	996798
Filling Instrument, Hollow Handle, Blue Teflon Coating, Type 5		2.6mm 2.2mm	996799
Filling Instrument, Hollow Handle, Blue Teflon Coating, Type 6		2.1mm 1.4mm	996800
Filling Instrument, Hollow Handle, Blue Teflon Coating, Type 7		1.1mm 1.3mm	996801
Filling Instrument, Hollow Handle, Blue Teflon Coating, Type 8		1.4mm 2.0mm 1.3mm	996802

INSTRUMENTS












INSTRUMENTS

NuMedical

Filling Instruments, Double-Ended, Teflon Coating

\$23.95 each

Code

Filling Instrument, Blue Teflon Coating, Type 1	 1.5mm	2.4mm 996803
Filling Instrument, Blue Teflon Coating, Type 2	 2.4mm	1.5mm 996804
Filling Instrument, Blue Teflon Coating, Type 3	 1.8mm	1.4mm 996805
Filling Instrument, Blue Teflon Coating, Type 4	 2.5mm	1.8mm 996806
Filling Instrument, Blue Teflon Coating, Type 5	 1.9mm	1.9mm 996807
Filling Instrument, Blue Teflon Coating, Type 6	 1.9mm 2.3mm	1.9mm 2.3mm 996808
Filling Instrument, Blue Teflon Coating, Type 7	 2.0mm	2.1mm 2.8mm 996809
Filling Instrument, Blue Teflon Coating, Type 8	 2.2mm 2.7mm	2.2mm 2.7mm 996810
Filling Instrument, Blue Teflon Coating, Type 9	 1.8mm	1.8mm 996811
Filling Instrument, Blue Teflon Coating, Type 10	 2.1mm 2.6mm	2.2mm 2.7mm 996812
Filling Instrument, Blue Teflon Coating, Type 11	 2.1mm 2.6mm	2.1mm 2.6mm 996813

Lining Applicator, Double-Ended

\$5.95 each

Code

Lining Applicator Double-Ended		996681
-----------------------------------	--	--------

Burnisher, 21B

\$4.95 each

Code

BURNISHER, 21B	 2.1mm 2.8mm	996842
----------------	---	--------

Bone File

\$16.95 each

Code

Bone File	 7.5mm 5.0mm 5.5mm	996874
-----------	--	--------

INSTRUMENTS

INSTRUMENTS

NuMedical

Amalgam Instruments ,Double-Ended,

\$5.95 each

Code

walls Carver,3



996776

Amalgam Instruments ,Double-Ended,

\$5.95 each

Code

Burnisher Ball,
Double-Ended,
27&29



996840

Explorer, Double Ended

\$10.95 each

Code

Stainless Steel Ex-
plorer, Double Ended



996542

Surgical Bone Curette

\$22.95 each

Code

Surgical Bone Curette



996875

Trimmer Mitchell, NO.4

\$4.95 each

Code

Trimmer Mitchell, No.4



996898

Cleoid Discoid Carver, Large

\$4.95 each

Code

Cleoid Discoid Carver
Large



996839

Periosteal Elevator Molt, NO.9

\$12.95 each

Code



996557

The rounded edge of the periosteal instrument minimises tissue tearing when lifting the gum flap while the sharper end makes the initial lifting of the flap easier.

INSTRUMENTS

INSTRUMENTS & ACCESSORIES

NuMedical



\$91.90
each

LED Light Mirror

Designed for clearer examination of the oral cavity. 2 AA batteries included. Mirror is removable and autoclavable up to 135°C.

Code	Colour
993231	Silver



\$355.50
/set

LED Light Mirror Set

Long life of 20,000 hours. Long handle to avoid infection with anti-slip silicon grip. Each set comes with 1 LED lamp, 10 mirrors and 5 AA batteries. Mirror is removable and autoclavable up to 135°C.

Code	Colour
993230	Silver

2, 4, 6, 8



\$18.85
/bag

10pcs/bag

Periodontal Probe and Explorer

This periodontal resin probe has a scale of 2mm, 4mm, 6mm and 8mm on one end. The resin probe is bendable and flexible. On another end of the probe is a #6 stainless steel explorer. Autoclavable up to 135°C.

Code	Colour
993269	Cream

3, 6, 9, 12



\$19.95
/bag

10pcs/bag

Periodontal Probe and Explorer

This periodontal resin probe has a scale of 3mm, 6mm, 9mm and 12mm on one end. The resin probe is bendable and flexible. On another end of the probe is a #6 stainless steel explorer. Autoclavable up to 135°C.

Code	Colour
993270	Cream



\$23.00
/box

50kits/box



Sterile Instrument Kit, 2pcs

ETO gas sterilised instrument kit containing 1 probe and 1 mirror per kit.

Code
990730



\$38.35
/box

50kits/box




Sterile Instrument Kit, 3pcs

ETO gas sterilised instrument kit containing 1 probe, 1 mirror and 1 tweezer per kit.

Code
990731

INSTRUMENTS

INSTRUMENTS




\$9.20
each

Amalgam Carrier - Straight

Ideal for carrying dental amalgam, a combination of liquid mercury and metal alloy mixture for filling tooth cavities. It has a push handle, angled deliver socket and spring bolt. Autoclavable up to 135°C. Reusable. Non-sterile.

Code	Colour
990881	Blue
990884	White



\$9.20
each

Amalgam Carrier - Mini Right

Ideal for carrying dental amalgam, a combination of liquid mercury and metal alloy mixture for filling tooth cavities. It has a push handle, angled deliver socket and spring bolt. Autoclavable up to 135°C. Reusable. Non-sterile.

Code	Colour
990882	Blue
990885	White



\$9.20
each

Amalgam Carrier - Right Angle

Ideal for carrying dental amalgam, a combination of liquid mercury and metal alloy mixture for filling tooth cavities. It has a push handle, angled deliver socket and spring bolt. Autoclavable up to 135°C. Reusable. Non-sterile.

Code	Colour
990883	Blue
990886	White



\$6.95
each

Molar Band Seater

Molar Band Seater. Autoclavable up to 135°C.

Code
996846



\$19.95
/box

20pcs/box

Autoclavable Mouth Mirrors

Double-sided mouth mirror with clear visibility. Autoclavable up to 135°C. Size: 175mm (L)

Code	Colour
993233	Blue
993234	Pink
993235	Green
993236	Yellow
993237	Purple



\$3.95
/Box

10pcs/box

Disposable Mouth Mirror

Disposable anti-fog fibre-glass mouth mirror. Mirror size: 22mm (Diameter). Length: 167mm.

Code
993245

INSTRUMENTS CLEAN

INSTRUMENTS CLEAN

NuMedical



\$99.95
each

Instrument Cassette, 20 Instruments

Stainless steel instrument cassette with closure and silicone inserts. Suitable for 20 instruments. Autoclavable up to 135°C. Size: 280mm(L) x 195mm(W) x 32mm(H).

Code
996820



\$59.95
each

Instrument Cassette, 10 Instruments

Stainless steel instrument cassette with closure and silicone inserts. Suitable for 10 instruments. Autoclavable up to 135°C. Size: 180mm(L) x 133mm(W) x 30mm(H).

Code
996819



\$39.95
each

Instrument Cassette, 7 Instruments

Stainless steel instrument cassette with closure and silicone inserts. Suitable for 7 instruments. Autoclavable up to 135°C. Size: 187mm(L) x 97mm(W) x 25mm(H).

Code
996818



\$33.95
each

Instrument Cassette, 5 Instruments

Stainless steel instrument cassette with closure and silicone inserts. Suitable for 5 instruments. Autoclavable up to 135°C. Size: 180mm(L) x 60mm(W) x 20mm(H).

Code
996817




\$26.95
each

Instrument Tray, with lid

Instrument Tray, with lid. Autoclavable up to 135°C. Size: 230mm(L) x 130mm(W) x 80mm(H).

Code
996822



\$4.95
each


Instrument Scaler Tray, Grooved Base

Instrument Scaler Tray, Grooved Base. Autoclavable up to 135°C. Size: 206mm(L) x 117mm(W) x 18mm(H).

Code
996821

INSTRUMENTS CLEAN

INSTRUMENTS CLEAN



\$199.00
each

Instrument Cassette, 16 Instruments with Space

Stainless steel instrument cassette with closure and silicone inserts. Suitable for 16 elevators. Autoclavable up to 135°C.
Size: 370mm(L) x 212mm(W) x 37mm(H).

Code
996838



\$41.95
each

Soaker Holder

Convenient: Use with cleaners, disinfectants or cold sterilants. Easy-to Use: Inner basket is easily removed and allows for convenient draining while hinged lid contains calendar to notate when solution needs to be changed.

High-Capacity: Holds up to 1 gallon of solution.
Size: 35.2cm(L) x 24.4cm(W) x 13cm(H).

Code
996883



\$16.45
/bag

Instrument Tie

Used to tie instruments together during sterilisation. Autoclavable up to 135°. Size: 100mm long.

Code	Colour
992968	Blue
992969	Red
992970	Green
992971	Yellow
992972	Purple

NuMedical



iSmile General Purpose Cleaner, Made in USA

Formulated to remove of pumice, rouge, tripoli, buffing compounds and other foreign matter from dentures, bridges, crowns, orthodontic appliances, instruments, diamonds, burs, etc. **Directions:** Shake before use. Mix 1 part General Purpose Cleaner with 10 parts water, immerse items directly into tank for 2-10 minutes. Rinse thoroughly. Makes 11 gallons. Must be used as the ultrasonic conductor when using chemicals in glass beakers. **Contents:** Non-ammoniated. Contains a rust inhibitor. Concentrated detergent mixture contains nonionic surfactants. **Cautions:** For Professional Use Only. Not for Intrnal Use. Keep Away From Children. **Warning:** Causes serious eye irritation. Harmful if Swallowed. Wear Protective Gloves, Protective Clothing and Eye Protection.

\$41.95
each

4L/bottle



View Details

Code
996879

INSTRUMENTS CLEAN

INSTRUMENTS CLEAN

NuMedical



iSmile Tartar/Stain Remover Ultrasonic Solution, Made in USA

Formulated for the removal of calculus, tartar, tobacco and other food stains from dentures, bridges, orthodontic appliances, etc. Also can be used for the removal of temporary and permanent cements. **Directions:** Pour this ready to use solution into a covered glass beaker and immerse in General Purpose Cleaner. The General Purpose Cleaner acts as the ultrasonic conductor between the beaker and the tank. Cover item in beaker with adequate amount of Tartar & Stain Remover and clean for 5-10 minutes. Items cleaned should be rinsed and neutralized in a solution of 1 easpoon bicarbonate of soda diluted to 1 cup of water.

Caution: Do not use on aluminum. Do not use directly in ultrasonic tank. For Professional Use Only. Not For Internal Use. Keep Away From Children. **Warning:** Causes severe skin burne and eye damage. Harmful if swallowed. Wear protective gloves, protective clothing and eye protection.

\$41.95
each

4L/bottle



View Details

Code
996880



iSmile Enzymatic Tablets, for Cleaning, Soaking, Ultrasonic, General Purpose, Evacuation Cleaner, Made in USA

Direction For Use: **1) For Cleaning Purpose:** Place two tablets for a typical 3.8Lt. It is preferable to use warm tap water not to exceed 60°C. Manually clean instruments with residual debris with a handled brush. Instruments should be sonicated for a minimun of 10 minutes. For tough stains add 5 or more minutes to the procedure. Rinse instruments with water to remove any excess product residue, dry and start your sterilization process. **2) For Soaking:** Place two tablets for a typical 3.8Lt. It is preferable to use warm tap water not to exceed 60°C. Soaking instruments overnight as a daily procedure is not recommended; However if this pcedure is preferred, the anti-corrosive formula will protect the instrument sfrom any corrosion. **3) For Ultrasonic:** Place two tablets for a typical 3.8Lt. It is preferable to use warm tap water not to exceed 60°C. Ultrasonic solution should be replaced daily. Do not add tables to any other disinfecting products. **4) General purpose Cleaner:** Place two tablets per 6-8Lt of water. It is preferable to use warm tap water not to exceed 60°C. For light cleaning around the office. **5) For Evacuation Cleaner:** Place two tablets per 3.8Lt. It is preferable to use warm tap water not to exceed 60°C. Can be used as a secondary evacuation cleaning agent. However it is recommended to use a shock solution/tablet as your primary evacuation cleaner due to its high concentration level and various effective cleaning agents that are not found in iSmile Ultrasonic Tablets. **6) Special:** Non-corrosive, non-ionic and biodegradable. Mint Fragrance. • 2 tablets yields 1 gallon(3.8Lt) of solution. **7) Warning:** Wear gloves and eye protection when hanling this product. In the event of contact with skin or eyes, wash immediately with water. Keep out of reach of children. Not for human consumption.

\$69.95
/box

64 pcs/box



View Details

Code
996881



iSmile Evacuation System Cleaner, up to 34 Cleanings, Made in USA

Triple Enzyme System; Non-Foaming Formula; Deodorizes as it cleans; Dissolves Organic Fluids and Debris; Environmentally Friendly, Neutral PH Formula; Safe for use with Amalgam Separators; Daily Use Improves Cleaning Efficiency; Mint-scent; **Mixing Directions:** Add 2 ounces(60 ml) of iSmile Evacuation System Cleaner to 30 ounces(900 ml) of water. Mix well and draw through evacuation system, high volume and saliva ejection hoses. **First Aid:** In case of contact, immediately flush eyesLitre or skin with plenty of water for at least 15 min. If solution is swallowed or comes in contact with eyes immediately seek immediate meical attention. **Caution:** Not for human consumption. Harmful if swallowed. Do not get in eyes, skin or clothing. For Professional Use Only. Keep Away From Children.

\$41.95
each

2L/bottle



View Details

Code
997855

INSTRUMENTS CLEAN



Professional Ultrasonic Cleaner

\$999.00
each

A professional ultrasonic cleaner! Features: **1. Large Tank Capacity:** Tank capacity is 7.5L(2.0gallon) and tank size is 33.1cm x 23.8cm x 10.5cm(13" x 9.4" x 4.1"). **2. Control Panel with Capacitive Sensing Technology:** Durable, Reliable, Resistant to water and harsh chemicals. **3. Drainage Valve:** Convenient to use. **4. Three Industrial Grade Ultrasonic Transducers:** Uniform distribution of ultrasonic waves, strong cleaning, high durability. **5. Three Ceramic Heaters:** Mounted underneath the bottom of the tank for safe and even heating and durable construction. **6. 3-Color LED Display:** Easy to operate, with 4 key control panel, 6 timer settings, 3 temperature settings. **7. Independent Control Circuit for Each Ultrasonic Transducer:** 40% higher efficiency, not susceptible to interference, more durable. **8. Multiple Circuit Protectors:** When overloaded or improperly used, the protectors shut down the power to certain areas to protect the machine. **9. Moisture-proofed PCB with Industrial IC:** PCBs are placed vertically for better moisture-proof performance. Better anti-interference. **10. Three Cooling Fans:** Two for cooling of PCB boards, one for cooling of ultrasonic transducers. Improve heat dissipation and beneficial for continuous operations. **11. Multiple Optional Accessories:** Suspensible plastic basket, stainless steel rack and tray, indirect cleaning tank, beaker and double beaker holder. **12. Solid and Advanced Structural Design:** 4.5mm thick wall housing made from engineering grade plastic with tongue in groove design for better water-proof and drop-proof performances. 60% quieter than steel housing designs, suitable for indoor uses. **(Accessories are not include in the cleaner)**



[View Details](#)

Code
997036

Code	Product		Price
997033	Stainless Steel Rack		\$29.90 each
997034	Indirect Cleaning Tank		\$29.90 each
997035	Stainless Steel Tray		\$69.90 each
997036	Beaker and Beaker Holder		\$89.90 each

INSTRUMENTS CLEAN



INSTRUMENTS CLEAN

Nu^omedical