

ORTHODONTICS

ORTHODONTICS



\$8.95
/tube

30pcs/tube

Stainless Steel Straight Wires Round

Size	Code
0.012"	997740
0.014"	997741
0.016"	997742
0.018"	997743
0.020"	997744



\$12.95
/tube

30pcs/tube

Australia Straight Wires Round

Size	Code
0.014"	997755
0.016"	997756
0.018"	997757
0.020"	997758




Stainless Steel Straight Wires Rectangular

\$16.95
/tube

30pcs/tube

Size	Code
0.016"x0.016"	997745
0.016"x0.022"	997746
0.014"x0.025"	997747
0.016"x0.025"	997748
0.017"x0.022"	997749
0.017"x0.025"	997750
0.018"x0.022"	997751
0.018"x0.025"	997752
0.019"x0.025"	997753
0.021"x0.025"	997754

NuMedical




Stainless Steel Archwires, Universal, Gold-Plated

10pcs/pack

Gold-Plated Stainless Steel archwires, universal form. Available in round and rectangle.

Size	Round	Price
0.012"	997067G	\$7.95/pack
0.014"	997068G	\$7.95/pack
0.016"	997069G	\$7.95/pack
0.018"	997070G	\$7.95/pack
0.020"	997071G	\$7.95/pack

Size	Rectangle	Price
0.016" x 0.016"	997072G	\$13.35/pack
0.016" x 0.022"	997073G	\$13.35/pack
0.017" x 0.022"	997074G	\$13.35/pack
0.017" x 0.025"	997075G	\$13.35/pack
0.018" x 0.018"	997076G	\$13.35/pack
0.018" x 0.022"	997077G	\$13.35/pack
0.018" x 0.025"	997078G	\$13.35/pack
0.019" x 0.025"	997079G	\$13.35/pack
0.020" x 0.025"	997080G	\$13.35/pack
0.021" x 0.025"	997081G	\$13.35/pack




Stainless Steel Archwires, Universal, Bright

10pcs/pack

Stainless Steel archwires with bright finish, universal form. Available in round and rectangle.

Size	Round	Price
0.012"	997067B	\$5.45/pack
0.014"	997068B	\$5.45/pack
0.016"	997069B	\$5.45/pack
0.018"	997070B	\$5.45/pack
0.020"	997071B	\$5.45/pack

Size	Rectangle	Price
0.016" x 0.016"	997072B	\$9.95/pack
0.016" x 0.022"	997073B	\$9.95/pack
0.017" x 0.022"	997074B	\$9.95/pack
0.017" x 0.025"	997075B	\$9.95/pack
0.018" x 0.018"	997076B	\$9.95/pack
0.018" x 0.022"	997077B	\$9.95/pack
0.018" x 0.025"	997078B	\$9.95/pack
0.019" x 0.025"	997079B	\$9.95/pack
0.020" x 0.025"	997080B	\$9.95/pack
0.021" x 0.025"	997081B	\$9.95/pack




Stainless Steel Archwires, Natural, Gold-Plated 10pcs/pack

Gold-Plated Stainless Steel archwires, natural form. Available in round and rectangle, and in Upper Right (UR) and Lower Left (LL).

Size	Round	Price
0.012", UR	997082G	\$7.95/pack
0.012", LL	997083G	\$7.95/pack
0.014", UR	997084G	\$7.95/pack
0.014", LL	997085G	\$7.95/pack
0.016", UR	997086G	\$7.95/pack
0.016", LL	997087G	\$7.95/pack
0.018", UR	997088G	\$7.95/pack
0.018", LL	997089G	\$7.95/pack
0.020", UR	997090G	\$7.95/pack
0.020", LL	997091G	\$7.95/pack

Size	Rectangle	Price
0.016" x 0.016", UR	997092G	\$13.35/pack
0.016" x 0.016", LL	997093G	\$13.35/pack
0.016" x 0.022", UR	997094G	\$13.35/pack
0.016" x 0.022", LL	997095G	\$13.35/pack
0.017" x 0.022", UR	997096G	\$13.35/pack
0.017" x 0.022", LL	997097G	\$13.35/pack
0.017" x 0.025", UR	997098G	\$13.35/pack
0.017" x 0.025", LL	997099G	\$13.35/pack
0.018" x 0.018", UR	997300G	\$13.35/pack
0.018" x 0.018", LL	997301G	\$13.35/pack
0.018" x 0.022", UR	997302G	\$13.35/pack
0.018" x 0.022", LL	997303G	\$13.35/pack
0.018" x 0.025", UR	997304G	\$13.35/pack
0.018" x 0.025", LL	997305G	\$13.35/pack
0.019" x 0.025", UR	997306G	\$13.35/pack
0.019" x 0.025", LL	997307G	\$13.35/pack
0.020" x 0.025", UR	997308G	\$13.35/pack
0.020" x 0.025", LL	997309G	\$13.35/pack
0.021" x 0.025", UR	997310G	\$13.35/pack
0.021" x 0.025", LL	997311G	\$13.35/pack




Stainless Steel Archwires, Natural, Bright 10pcs/pack

Stainless Steel archwires with bright finish, natural form. Available in round and rectangle, and in Upper Right (UR) and Lower Left (LL).

Size	Round	Price
0.012", UR	997082B	\$5.45/pack
0.012", LL	997083B	\$5.45/pack
0.014", UR	997084B	\$5.45/pack
0.014", LL	997085B	\$5.45/pack
0.016", UR	997086B	\$5.45/pack
0.016", LL	997087B	\$5.45/pack
0.018", UR	997088B	\$5.45/pack
0.018", LL	997089B	\$5.45/pack
0.020", UR	997090B	\$5.45/pack
0.020", LL	997091B	\$5.45/pack

Size	Rectangle	Price
0.016" x 0.016", UR	997092B	\$9.95/pack
0.016" x 0.016", LL	997093B	\$9.95/pack
0.016" x 0.022", UR	997094B	\$9.95/pack
0.016" x 0.022", LL	997095B	\$9.95/pack
0.017" x 0.022", UR	997096B	\$9.95/pack
0.017" x 0.022", LL	997097B	\$9.95/pack
0.017" x 0.025", UR	997098B	\$9.95/pack
0.017" x 0.025", LL	997099B	\$9.95/pack
0.018" x 0.018", UR	997300B	\$9.95/pack
0.018" x 0.018", LL	997301B	\$9.95/pack
0.018" x 0.022", UR	997302B	\$9.95/pack
0.018" x 0.022", LL	997303B	\$9.95/pack
0.018" x 0.025", UR	997304B	\$9.95/pack
0.018" x 0.025", LL	997305B	\$9.95/pack
0.019" x 0.025", UR	997306B	\$9.95/pack
0.019" x 0.025", LL	997307B	\$9.95/pack
0.020" x 0.025", UR	997308B	\$9.95/pack
0.020" x 0.025", LL	997309B	\$9.95/pack
0.021" x 0.025", UR	997310B	\$9.95/pack
0.021" x 0.025", LL	997311B	\$9.95/pack




Stainless Steel Archwires, Standard, Gold-Plated

10pcs/pack

Gold-Plated Stainless Steel archwires, standard form. Available in round and rectangle, and in Upper Right (UR) and Lower Left (LL).

Size	Round	Price
0.012", UR	997312G	\$7.95/pack
0.012", LL	997313G	\$7.95/pack
0.014", UR	997314G	\$7.95/pack
0.014", LL	997315G	\$7.95/pack
0.016", UR	997316G	\$7.95/pack
0.016", LL	997317G	\$7.95/pack
0.018", UR	997318G	\$7.95/pack
0.018", LL	997319G	\$7.95/pack
0.020", UR	997320G	\$7.95/pack
0.020", LL	997321G	\$7.95/pack

Size	Rectangle	Price
0.016" x 0.016", UR	997322G	\$13.35/pack
0.016" x 0.016", LL	997323G	\$13.35/pack
0.016" x 0.022", UR	997324G	\$13.35/pack
0.016" x 0.022", LL	997325G	\$13.35/pack
0.017" x 0.022", UR	997326G	\$13.35/pack
0.017" x 0.022", LL	997327G	\$13.35/pack
0.017" x 0.025", UR	997328G	\$13.35/pack
0.017" x 0.025", LL	997329G	\$13.35/pack
0.018" x 0.018", UR	997330G	\$13.35/pack
0.018" x 0.018", LL	997331G	\$13.35/pack
0.018" x 0.022", UR	997332G	\$13.35/pack
0.018" x 0.022", LL	997333G	\$13.35/pack
0.018" x 0.025", UR	997334G	\$13.35/pack
0.018" x 0.025", LL	997335G	\$13.35/pack
0.019" x 0.025", UR	997336G	\$13.35/pack
0.019" x 0.025", LL	997337G	\$13.35/pack
0.020" x 0.025", UR	997338G	\$13.35/pack
0.020" x 0.025", LL	997339G	\$13.35/pack
0.021" x 0.025", UR	997340G	\$13.35/pack
0.021" x 0.025", LL	997341G	\$13.35/pack



Stainless Steel Archwires, Standard, Bright

10pcs/pack

Stainless Steel archwires with bright finish, standard form. Available in round and rectangle, and in Upper Right (UR) and Lower Left (LL).

Size	Round	Price
0.012", UR	997312B	\$5.45/pack
0.012", LL	997313B	\$5.45/pack
0.014", UR	997314B	\$5.45/pack
0.014", LL	997315B	\$5.45/pack
0.016", UR	997316B	\$5.45/pack
0.016", LL	997317B	\$5.45/pack
0.018", UR	997318B	\$5.45/pack
0.018", LL	997319B	\$5.45/pack
0.020", UR	997320B	\$5.45/pack
0.020", LL	997321B	\$5.45/pack

Size	Rectangle	Price
0.016" x 0.016", UR	997322B	\$9.95/pack
0.016" x 0.016", LL	997323B	\$9.95/pack
0.016" x 0.022", UR	997324B	\$9.95/pack
0.016" x 0.022", LL	997325B	\$9.95/pack
0.017" x 0.022", UR	997326B	\$9.95/pack
0.017" x 0.022", LL	997327B	\$9.95/pack
0.017" x 0.025", UR	997328B	\$9.95/pack
0.017" x 0.025", LL	997329B	\$9.95/pack
0.018" x 0.018", UR	997330B	\$9.95/pack
0.018" x 0.018", LL	997331B	\$9.95/pack
0.018" x 0.022", UR	997332B	\$9.95/pack
0.018" x 0.022", LL	997333B	\$9.95/pack
0.018" x 0.025", UR	997334B	\$9.95/pack
0.018" x 0.025", LL	997335B	\$9.95/pack
0.019" x 0.025", UR	997336B	\$9.95/pack
0.019" x 0.025", LL	997337B	\$9.95/pack
0.020" x 0.025", UR	997338B	\$9.95/pack
0.020" x 0.025", LL	997339B	\$9.95/pack
0.021" x 0.025", UR	997340B	\$9.95/pack
0.021" x 0.025", LL	997341B	\$9.95/pack

ORTHODONTICS

ORTHODONTICS



Super Elastic NiTi Archwires - Round

Super elastic round archwires made from Nickel Titanium. Available in square, ovoid and nature for upper and lower arch.

\$4.45 /pack

10pcs/pack

Size	UPPER			LOWER		
	Square	Ovoid	Nature	Square	Ovoid	Nature
0.012"	995942	995952	995962	995947	995957	995967
0.014"	995943	995953	995963	995948	995958	995968
0.016"	995944	995954	995964	995949	995959	995969
0.018"	995945	995955	995965	995950	995960	995970
0.020"	995946	995956	995966	995951	995961	995971



Super Elastic NiTi Archwires - Rectangular

Super elastic rectangular archwires made from Nickel Titanium. Available in square, ovoid and nature for upper and lower arch.

\$7.55 /pack

10pcs/pack

Size	UPPER			LOWER		
	Square	Ovoid	Nature	Square	Ovoid	Nature
0.016" x 0.016"	995972	995988	996004	995980	995996	996012
0.016" x 0.022"	995973	995989	996005	995981	995997	996013
0.017" x 0.022"	995974	995990	996006	995982	995998	996014
0.017" x 0.025"	995975	995991	996007	995983	995999	996015
0.018" x 0.022"	995976	995992	996008	995984	996000	996016
0.018" x 0.025"	995977	995993	996009	995985	996001	996017
0.019" x 0.025"	995978	995994	996010	995986	996002	996018
0.021" x 0.025"	995979	995995	996011	995987	996003	996019

NuMedical

ORTHODONTICS



Heat Activated NiTi Archwires - Round

Heat-activated round archwires made from Nickel Titanium.
Available in square, ovoid and nature for upper and lower arch.

\$3.95
/pack

5pcs/pack



Size	UPPER			LOWER		
	Square	Ovoid	Nature	Square	Ovoid	Nature
0.012"	996020	996030	996040	996025	996035	996045
0.014"	996021	996031	996041	996026	996036	996046
0.016"	996022	996032	996042	996027	996037	996047
0.018"	996023	996033	996043	996028	996038	996048
0.020"	996024	996034	996044	996029	996039	996049



Heat Activated NiTi Archwires - Rectangular

Heat-activated rectangular archwires made from Nickel Titanium.
Available in square, ovoid and nature for upper and lower arch.

\$5.95
/pack

5pcs/pack



Size	UPPER			LOWER		
	Square	Ovoid	Nature	Square	Ovoid	Nature
0.016" x 0.016"	996050	996066	996082	996058	996074	996090
0.016" x 0.022"	996051	996067	996083	996059	996075	996091
0.017" x 0.022"	996052	996068	996084	996060	996076	996092
0.017" x 0.025"	996053	996069	996085	996061	996077	996093
0.018" x 0.022"	996054	996070	996086	996062	996078	996094
0.018" x 0.025"	996055	996071	996087	996063	996079	996095
0.019" x 0.025"	996056	996072	996088	996064	996080	996096
0.021" x 0.025"	996057	996073	996089	996065	996081	996097

ORTHODONTICS

NuMedical

ORTHODONTICS

ORTHODONTICS

Nu2medical



\$3.95
/box

10 pcs/box

UDG Stainless Steel Archwires Ovoid Rectangular

Size	Upper	Lower
0.016" x 0.016"	996225	996226
0.016" x 0.022"	996227	996228
0.014" x 0.025"	996229	996230
0.016" x 0.025"	996231	996232
0.017" x 0.022"	996233	996234
0.017" x 0.025"	996235	996236
0.018" x 0.022"	996237	996238
0.018" x 0.025"	996239	996240
0.019" x 0.025"	996241	996242
0.021" x 0.025"	996243	996244



\$3.95
/box

10 pcs/box

UDG Stainless Steel Archwires Ovoid Round

Size	Upper	Lower
0.012"	996215	996216
0.014"	996217	996218
0.016"	996219	996220
0.018"	996221	996222
0.020"	996223	996224



\$3.95
/box

10 pcs/box

UDG NiTi Super Elastic Arch- wires Natural Rectangular

Size	Upper	Lower
0.016" x 0.016"	996255	996256
0.016" x 0.022"	996257	996258
0.014" x 0.025"	996259	996260
0.016" x 0.025"	996261	996262
0.017" x 0.022"	996263	996264
0.017" x 0.025"	996265	996266
0.018" x 0.022"	996267	996268
0.018" x 0.025"	996269	996270
0.019" x 0.025"	996271	996272
0.021" x 0.025"	996273	996274



\$3.95
/box

10 pcs/box

UDG NiTi Super Elastic Arch- wires Natural Round

Size	Upper	Lower
0.012"	996245	996246
0.014"	996247	996248
0.016"	996249	996250
0.018"	996251	996252
0.020"	996253	996254

ORTHODONTICS



Copper Niti Archwires Ovoid Round

\$32.95
/box

10 pcs/box

Size	Upper	Lower
0.012"	997710	997711
0.014"	997712	997713
0.016"	997714	997715
0.018"	997716	997717
0.020"	997718	997719



Copper Niti Archwires Ovoid Rectangular

\$32.95
/box

10 pcs/box

Size	Upper	Lower
0.016" x0.016"	997720	997721
0.016" x0.022"	997722	997723
0.014" x0.025"	997724	997725
0.016" x0.025"	997726	997727
0.017" x0.022"	997728	997729
0.017" x0.025"	997730	997731
0.018" x0.022"	997732	997733
0.018" x0.025"	997734	997735
0.019" x0.025"	997736	997737
0.021" x0.025"	997738	997739

ORTHODONTICS

Nu
medical



IPR STRIP

1. Ergonomic, single-handed design provides optimal tactile control and reduces hand fatigue.
2. Eliminates the need to hold small restorations and use rotary instruments or articulating films.
3. Central opening for enhanced visual perception and access for tools.
4. Will not remove excess enamel or cut lips or gums.
5. Flexible strips conform to the natural contours of the teeth to avoid creating sharp corners.
6. Autoclavable and reusable.

\$78.95
/box

10 pcs/box

Code		Color	Size
997796		Black	strip width 3mm, thickness 0.1mm
997797		Yellow	strip width 4mm, thickness 0.08mm, double-sided & ultra-fine diamond
997798		Dark Blue	strip width 4mm, thickness 0.1mm, double-sided & ultra-fine diamond
997799		Light Blue	strip width 4mm, thickness 0.15mm, single-sided & extra-fine diamond
997800		Purple	strip width 4mm, thickness 0.2mm, single-sided & medium-coarse diamond
997801		Red	strip width 4mm, thickness 0.2mm, double-sided & extra-fine diamond
997802		White	strip width 4mm, thickness 0.3mm, double-sided & coarse diamond
997803		Orange	strip width 4mm, thickness 0.3mm, double-sided & medium coarse diamond,
997804		Pink	strip width 4mm, thickness 0.4mm, double-sided & coarse diamond,
997805		Green	strip width 4mm, thickness 0.01mm, single-sided & extra-fine diamond,

ORTHODONTICS

ORTHODONTICS

NuMedical



\$29.95
/box
20 pcs/box



UDG Ceramic Brackets 4/5 with Hooks, 0.22

Size	Code
ROTH	996275
MBT	996278



\$8.95
/box
20 pcs/box



UDG Metal Brackets 3/4/5 with Hooks Mini

Size	Code
ROTH	996209
MBT	996212



\$199.95
/box
20 pcs/box

Sapphire Bracket with Position Cross Roth 0.022 345hook

Lowest friction aesthetic sapphire ceramic bracket. Superior crystal clear appearance. Low profile design for maximum patient comfort. Mesh base design ensure better bonding strength.

Size	Code
ROTH	996151
MBT	996152



\$359.00
/box
20pcs/box

UDG Ceramic Active Self-Ligating Brackets Mini 45 with Hooks 0.22

Size	Code
ROTH	996276
MBT	996279



\$75.50
/box
20pcs/box

UDG Metal Passive Self-Ligating Brackets, 0.22

Size	Code
ROTH Mini 345 with Hooks	996211



\$359.00
/box
20pcs/box

UDG Ceramic Passive Self-Ligating Brackets Mini 45 with Hooks 0.22

Size	Code
ROTH	996277
MBT	996280

ORTHODONTICS

ORTHODONTICS

Nu2medical



Mini Ceramic Brackets

- Pure Polycrystalline aesthetic ceramic
- Twin wing
- Contoured design with rounded edge for maximum patient comfort
- Optimal base for superior bonding/debonding
- Low profile height design

\$45.95
/set

20pcs/set

ROTH TYPE			MBT TYPE	
Size	With Hook 3	With Hook 3, 4, 5	With Hook 3	With Hook 3, 4, 5
0.018"	993362	993363	993366	993367
0.022"	993360	993361	993364	993365



Ultra Clear Mini Ceramic Brackets

- True twin wings and low profile design for excellent performance
- World's first 2-Step groove for easy and accurate positioning
- Rounded edge and full-body coating
- Plasma-silica coating for superior sliding and to minimise fracture risks
- Half base coating for superior retention and easy debonding
- Low profile height design

\$179.00
/set

20pcs/set

ROTH TYPE			MBT TYPE	
Size	With Hook 3	With Hook 3, 4, 5	With Hook 3	With Hook 3, 4, 5
0.018"	993454	993455	993458	993459
0.022"	993452	993453	993456	993457



Mini Ceramic Self-Ligating Brackets

- Rounded edge for maximum patient comfort
- Tooth-coloured clips with Rhodium coating
- Optimal base for superior bonding/debonding
- Low profile height design

\$359.00
/set

20pcs/set

ROTH TYPE	
Size	With Hook 3, 4, 5
0.022"	994395



\$19.95
each

Self-Ligating Bracket Tool

Assists with placement of brackets.

Code
994394



\$15.95
/pack

NiTi Open Springs

Nickel-titanium open coil springs deliver continuous force over extended ranges to open spaces to accommodate new teeth.

2 pcs/pack

	0.010"	0.020"
180mm	995880	995881



Mini Metal Brackets (Laser-Welded)

- Rounded edge for maximum patient comfort
- Low profile height design
- Optimal base for superior bonding/debonding
- Laser-welding technology

\$17.95
/set

20pcs/set

ROTH TYPE		
Size	With Hook 3	With Hook 3, 4, 5
0.018"	994392	994393
0.022"	994390	994391



Mini Metal Self-Ligating Brackets

- Rounded edge for maximum patient comfort
- Low profile height design
- Optimal base for superior bonding/debonding
- Self-ligating clip
- MIM advanced technology

\$145.00
/set

20pcs/set

ROTH TYPE		
Size	With Hook 3	With Hook 3, 4, 5
0.018"	994388	994389
0.022"	994386	994387



Standard Metal Brackets (80 Mesh Base)

- Rounded edge for maximum patient comfort
- Low profile height design
- Optimal base for superior bonding/debonding

\$27.95
/set

20pcs/set

ROTH TYPE				MBT TYPE		
Size	Without Hook	With Hook 3	With Hook 3, 4, 5	Without Hook	With Hook 3	With Hook 3, 4, 5
0.018"	994422	994423	994424	994428	994429	994430
0.022"	994419	994420	994421	994425	994426	994427



Mini Metal Brackets (80 Mesh Base)

- Rounded edge for maximum patient comfort
- Low profile height design
- Optimal base for superior bonding/debonding

\$27.95
/set

20pcs/set

ROTH TYPE				MBT TYPE		
Size	Without Hook	With Hook 3	With Hook 3, 4, 5	Without Hook	With Hook 3	With Hook 3, 4, 5
0.018"	994434	994435	994436	994440	994441	994442
0.022"	994431	994432	994433	994437	994438	994439

ORTHODONTICS



Buccal Tubes (Convertible)

- Low profile design reduces occlusal interference
- * Comfortable hooks for maximum patient comfort
- Anatomically conforms to teeth

\$1.65 each

MAXILLIARY SINGLE BONDABLE (RECTANGULAR TUBE)					ROTH TYPE		MBT TYPE	
	Torque	Distal Offset	Length	Designation	0.018"	0.022"	0.018"	0.022"
1st Molar	(-) 14°	10°	4	Upper Right	996193	996197	996195	996199
	(-) 14°	10°	4	Upper Left	996194	996198	996196	996200
MANDIBULAR SINGLE BONDABLE (RECTANGULAR TUBE)					MBT TYPE			
	Torque	Distal Offset	Length	Designation	0.018"	0.022"		
1st Molar	(-) 20°	0°	4	Lower Right	996203	996207		
	(-) 20°	0°	4	Lower Left	996204	996208		
MANDIBULAR SINGLE BONDABLE (RECTANGULAR TUBE)					ROTH TYPE			
	Torque	Distal Offset	Length	Designation	0.018"	0.022"		
1st Molar	(-) 25°	4°	4	Lower Right	996201	996205		
	(-) 25°	4°	4	Lower Left	996202	996206		



Buccal Tubes Trumpeted Opening (Non-Convertible)

- Low profile design enhances bonding strength
- * Smooth-surface hooks for maximum patient comfort
- Funneled entrance guides wire into the buccal tube

\$1.65 each

MAXILLIARY SINGLE BONDABLE (RECTANGULAR TUBE)					ROTH TYPE		MBT TYPE	
	Torque	Distal Offset	Length	Designation	0.018"	0.022"	0.018"	0.022"
1st Molar	(-) 14°	10°	4	Upper Right	995830	995834	995832	995836
	(-) 14°	10°	4	Upper Left	995831	995835	995833	995837
2nd Molar	(-) 14°	10°	4	Upper Right	995846	995850	995848	995852
	(-) 14°	10°	4	Upper Left	995847	995851	995849	995853
MANDIBULAR SINGLE BONDABLE (RECTANGULAR TUBE)					ROTH TYPE			
	Torque	Distal Offset	Length	Designation	0.018"	0.022"		
1st Molar	(-) 25°	4°	4	Lower Right	995838	995842		
	(-) 25°	4°	4	Lower Left	995839	995843		
2nd Molar	(-) 25°	4°	4	Lower Right	995854	996157		
	(-) 25°	4°	4	Lower Left	995855	996158		
MANDIBULAR SINGLE BONDABLE (RECTANGULAR TUBE)					MBT TYPE			
	Torque	Distal Offset	Length	Designation	0.018"	0.022"		
1st Molar	(-) 20°	0°	4	Lower Right	995840	995844		
	(-) 20°	0°	4	Lower Left	995841	995845		
2nd Molar	(-) 20°	0°	4	Lower Right	995856	996159		
	(-) 20°	0°	4	Lower Left	996156	996160		

ORTHODONTICS



Buccal Tubes (Non-Convertible)

- Low profile design reduces occlusal interference
- Comfortable hooks for maximum patient comfort
- Anatomically conforms to teeth

\$1.65
each

MAXILLIARY SINGLE BONDABLE (RECTANGULAR TUBE)					ROTH TYPE		MBT TYPE	
	Torque	Distal Offset	Length	Designation	0.018"	0.022"	0.018"	0.022"
1st Molar	(-) 14°	10°	4	Upper Right	996161	996165	996163	996167
	(-) 14°	10°	4	Upper Left	996162	996166	996164	996168
2nd Molar	(-) 14°	10°	4	Upper Right	996177	996181	996179	996183
	(-) 14°	10°	4	Upper Left	996178	996182	996180	996184

MANDIBULAR SINGLE BONDABLE (RECTANGULAR TUBE)					ROTH TYPE	
	Torque	Distal Offset	Length	Designation	0.018"	0.022"
1st Molar	(-) 25°	4°	4	Lower Right	996169	996173
	(-) 25°	4°	4	Lower Left	996170	996174
2nd Molar	(-) 25°	4°	4	Lower Right	996185	996189
	(-) 25°	4°	4	Lower Left	996186	996190

MANDIBULAR SINGLE BONDABLE (RECTANGULAR TUBE)					MBT TYPE	
	Torque	Distal Offset	Length	Designation	0.018"	0.022"
1st Molar	(-) 20°	0°	4	Lower Right	996171	996175
	(-) 20°	0°	4	Lower Left	996172	996176
2nd Molar	(-) 20°	0°	4	Lower Right	996187	996191
	(-) 20°	0°	4	Lower Left	996188	996192




Molar Bands, Pre-welded, Double Rectangular

- Pre-welded molar bands (0.022" Roth type)
- With double rectangular tube, convertible
- Each bag comes with 5pcs of Upper Right and 5pcs of Upper Left
- Available in 10 different sizes

\$25.95
/bag

10 pcs/bag

1ST MOLAR BAND SIZE (0.022" ROTH TYPE)									
36	36+	37	37+	38	38+	39	39+	40	40+
997342	997343	997344	997345	997346	997347	997348	997349	997350	997351



Molar Bands, Pre-welded, Double Convertible

- Pre-welded molar bands (0.022" Roth type)
- Double, convertible with head gear tube
- Each bag comes with 5pcs of Upper Right and 5pcs of Upper Left
- Available in 10 different sizes

\$25.95
/bag

10 pcs/bag

1ST MOLAR BAND SIZE (0.022" ROTH TYPE)									
36	36+	37	37+	38	38+	39	39+	40	40+
997352	997353	997354	997355	997356	997357	997358	997359	997360	997361

ORTHODONTICS

NuMedical

ORTHODONTICS

ORTHODONTICS



Molar Bands, Pre-welded, Single Convertible

- Pre-welded molar bands (0.022" Roth type)
- Single, convertible
- Each bag comes with 5pcs of Upper Right and 5pcs of Upper Left
- Available in 10 different sizes

\$25.95
/bag

10 pcs/bag

1ST MOLAR BAND SIZE (0.022" ROTH TYPE)									
36	36+	37	37+	38	38+	39	39+	40	40+
997362	997363	997364	997365	997366	997367	997368	997369	997370	997371



Molar Bands, Pre-welded, Single Non-Convertible

- Pre-welded molar bands (0.022" Roth type)
- Single, non-convertible
- Each bag comes with 5pcs of Upper Right and 5pcs of Upper Left
- Available in 10 different sizes

\$25.95
/bag

10 pcs/bag

1ST MOLAR BAND SIZE (0.022" ROTH TYPE)									
36	36+	37	37+	38	38+	39	39+	40	40+
997372	997373	997374	997375	997376	997377	997378	997379	997380	997381

NuMedical



\$36.95
each

2 pcs/bag

Double Ended Positioning Gauge

This gauge is placed on the incisal/occlusal surface of the tooth to simplify bracket positioning. Measures 2.0mm, 2.5mm, 3.0mm, 3.5mm and, 3.5mm, 4.0mm, 4.5mm 5.0mm.

Code
994476



\$13.75
each

Cross Bracket Positioning Gauge

This gauge is placed on the incisal/occlusal surface of the tooth to simplify bracket positioning. Measures 3.5mm, 4mm, 4.5mm and 5mm.

Code
996147



Heavy I



Heavy II



Light II Extra Long



Light II

\$17.55
pack

2 pcs/pack

J Hooks

Made from high quality stainless steel. Designed to pull anterior teeth inwards and backwards, and to help to move posterior teeth distally.

Width	Code
Heavy I, 85mm	995882
Heavy II, 85mm	995883
Light II Extra Long, 78mm	995884
Light II, 65mm	995885

ORTHODONTICS

ORTHODONTICS



Ball Top



T Top



Long Curved

Crimpable Hooks

10 pcs/pack

Made from high quality stainless steel. Easily place onto any archwire inside or outside the mouth. Secure anchor for any auxiliary attachment with bi-directional hook for easy ligation. Fits archwires from 0.022" (0.56mm) to 0.025" (0.64mm).

Width	Price	T Top	Ball Top	Long Curved
2.1mm	\$8.75/pack	995857	995862	
2.7mm	\$8.75/pack	995858	995863	
4.0mm	\$8.75/pack	995859	995864	
5.5mm	\$8.75/pack	995860	995865	
7.0mm	\$8.75/pack	995861	995866	
7.0mm (Universal left)	\$12.85/pack			995867
7.0mm (Universal right)	\$12.85/pack			995868



\$5.45 /pack

10 pcs/pack

Crimpable Stops

Crimpable stops can easily be placed anywhere on the archwire at a desired point. Will firmly remain in place during treatment
Width: Universal.

Code
995869



\$8.55 /pack

10 pcs/pack

NiTi Close Springs

Nickel-titanium closed coil springs deliver continuous force over extended ranges of activation and working times to assist maintaining space between two teeth.

Force	9mm	10mm	12mm	14mm	16mm
Light (0.008)	995870	995871	995872	995873	995874
Medium (0.01)	995875	995876	995877	995878	995879

NuMedical

Low Friction: It can reduce 70% friction coefficient between bracket and archwires through using new Teflon coating and metal surface polishing technology. The smaller the friction, the greater the force on the teeth and the shorter the correction time.

Quick Operation Design: The design of Self-Ligating and Active hooks supply quick operation in treatment,

Saving Time & Money: The reduction of friction coefficient will help professional dentists and patients to save almost 25% time of treatment.

Each box contains 12pcs Archwires.

NiTi Thermal Active, Round, 0.018" 2pcs(upper and lower), 0.016" 2pcs(upper and lower), 0.014" 2pcs(upper and lower).

NiTi Super Elastic, Rectangular, 0.017"x0.025" 2pcs(upper and lower), 0.016"x0.022" 2pcs(upper and lower).

Stainless Steel Archwires, Rectangular, 0.019"x0.025" 2pcs(upper and lower).

\$99.95 /pack



View Details



UDG Uni-slide & Super Low Friction Metal Passive Self-ligating Brackets With Hook 345 0.22" 20 Brackets + 12 Archwires

Code	Model	Size
997791	High Torqu	0.22"
997790	MBT	0.22"
997789	ROTH	0.22"
997792	Medium Torque	0.22"
997793	Low Torque	0.22"

ORTHODONTICS

ORTHODONTICS



Without hook and indent

Code
995905



Without hook with indent

Code
995906



With hook and indent
Size: Distal L

Code
995907



With hook and indent
Size: Distal R


Code
995908

\$12.30
/pack

10 pcs/pack

Lingual Sheaths

Stainless steel lingual sheaths for attaching accessories such as transpalatal arches NiTi molar rotators and expanders.



Short

Long

\$27.70
/pack

100 pcs/pack

Kobayashi Hooks

The long Kobayashi hook is used to place Class II or vertical elastics onto brackets, while the short Kobayashi pre-twisted hook is handy to work on rotated teeth.

Size	Code
Short	995887
Long	995886



Short Arm



Long Arm

\$6.95
/pack

10 pcs/pack

Lingual Cleats

A versatile cleat formed by stainless steel. Ends may be raised or lowered as required to engage auxiliaries. Edges and corners rounded for patient comfort.

Size	Code
Short Arm	995898
Long Arm	995899



\$6.55
/pack

2 pcs/pack

Lingual Retainers - Bondable

Bondable lingual splints which extend from cuspid to cuspid or can be shortened to retain any combination of adjacent teeth between cuspids.

Size	Lower	Size	Upper
29mm	995893	38mm	995888
31mm	995894	40mm	995889
33mm	995895	42mm	995890
35mm	995896	45mm	995891
37mm	995897	47mm	995892



Weldable



Bondable
Eyelets

\$12.00
/pack

10 pcs/pack



Bondable
Mesh Single



Bondable



Bondable
Mesh Wings

Lingual Buttons

To assist springs, wires and elastic ligatures for traction and teeth movement. Buttons can be fixed directly onto teeth using resins.

Size	Code	Type
4.0mm x 2.6mm	995900	Weldable
Ø 3.5mm	995901	Bondable Mesh Wings
3.48mm x 2.5mm	995902	Bondable Mesh Single
Ø 3.5mm	995903	Bondable Eyelets
Ø 3.0mm	995904	Bondable

NuMedical

\$3.45
/bag

4 pcs/bag



1st & 2nd Molar Bands

Molar Bands, 4pcs per bag, upper-right x 1pc, upper-left x 1pc, Lower-right x 1pc, Lower-left x 1pc. In our size list, 34 means 34mm perimeter of bands and 34+ means 34.5mm perimeter.

Size	Code	Size	Code
1st, 34	994443	2nd, 31	994458
1st, 34+	994444	2nd, 31+	994459
1st, 35	994445	2nd, 32	994460
1st, 35+	994446	2nd, 32+	994461
1st, 36	994447	2nd, 33	994462
1st, 36+	994448	2nd, 33+	994463
1st, 37	994449	2nd, 34	994464
1st, 37+	994450	2nd, 34+	994465
1st, 38	994451	2nd, 35	994466
1st, 38+	994452	2nd, 35+	994467
1st, 39	994453	2nd, 36	994468
1st, 39+	994454	2nd, 36+	994469
1st, 40	994455	2nd, 37	994470
1st, 40+	994456	2nd, 37+	994471
1st, 41	994457	2nd, 38	994472

\$87.95
each




Dry Field System

Eliminates saliva and fully exposes the patient's top and bottom arches making it ideal for all bonding applications.

Size	Colour	Code
Small	Red	994474
Large	White	994475

\$6.10
/roll

50gm/roll



Ligature Wire

50gm roll of ligature steel wire fastens archwires into brackets.

Diameter	Code
0.2mm	994477
0.25mm	994478
0.30mm	994479

\$54.95
each



Ortho-Model Storage Box, Top-Load

A practical and space-saving way to store orthodontics models. Each storage box holds up to 80 models. It has 4 tray layers. Each tray holds 20 models. Easy and quick to assemble.

Size: 515mm (L) x 495mm (W) x 316mm (H)

Code	Colour
993437	White

\$89.95
each



Ortho-Model Drawer System

A practical and space-saving way to store orthodontics models. Each storage box holds up to 80 models. It has 4 tray layers. Each tray holds 20 models. Easy and quick to assemble.

Size: 490mm (L) x 468mm (W) x 400mm (H)

Code	Colour
993438	White




Adapter Rings

\$139.65
each

Shoot XT-103C Macro LED Ring Flash

This Macro LED Ring Flash Light allows dental professionals to take clearer intraoral photos than ever before. With bright, HD-friendly soft light and flicker-free LED light output, unwanted shadows are eliminated. Fits camera lens with diameter of 52mm, 55mm, 58mm, 62mm, 67mm, 72mm and 77mm. Comes with 8pcs of adapter rings.

Code
996145




\$8.15
/spool

Power Chain - Continuous

4.572 metres/spool

Used to assist with teeth movement and commonly used together with brackets and archwires for extra force.

Code	Colour
993374	Clear
993375	White
993376	Blue
993377	Red
993378	Pink
993379	Yellow
993392	Lake Blue
993393	Silver
993394	Purple
993395	Green




\$8.15
/spool

Power Chain - Short

4.572 metres/spool

Used to assist with teeth movement and commonly used together with brackets and archwires for extra force.

Code	Colour
993380	Clear
993381	White
993382	Blue
993383	Red
993384	Pink
993385	Yellow
993396	Lake Blue
993397	Silver
993398	Purple
993399	Green



\$8.15
/spool

Power Chain - Long

4.572 metres/spool

Used to assist with teeth movement and commonly used together with brackets and archwires for extra force.

Code	Colour
993386	Clear
993387	White
993388	Blue
993389	Red
993390	Pink
993391	Yellow
993440	Lake Blue
993441	Silver
993442	Purple
993443	Green



\$8.15
/bag

40pcs/bag


Elastic Ligatures

Used to affix the archwire to brackets. Each stick comes with 26 ligatures. Each bag comes with 40 sticks.

Code	Colour
993403	Clear
993404	White
993405	Blue
993406	Red
993407	Pink
993408	Yellow
994059	Lake Blue
994060	Silver
994061	Purple
994062	Green
994063	Assorted

ORTHODONTICS

ORTHODONTICS




\$11.95
/bag

10pcs/bag

S-Shape Elastic Separator

Used to separate teeth, prior to placing molar bands. Each S-Shape stick comes with 70 individual separators. Each bag comes with 10 S-Shape sticks.

Code	Colour
996150	Blue



\$11.95
/bag

100pcs/bag


Elastic Ligature Directors

To easily remove and place elastic ligatures (colour ties) onto brackets. Assorted colours of green, yellow, red and blue.

Code	Colour
993373	Assorted

\$0.77
/bag

100pcs/bag



Elastics

Used to connect the top and bottom brackets with hooks. Elastics assist in adjusting the bite and jaw position. Force: Medium-Light (99.22 g).

Size	Code
1/8"	993368
3/16"	993369
1/4"	993370
5/16"	993371
3/8"	993372



\$37.95
each

Orthodontic Face Mask, Blue

Made from high quality stainless steel and plastic with leather padded contact points for maximum patient comfort. The forehead, chin pads and "cross piece" are adjustable. This device is used by patients with "underbite", and assists in realigning the patient's jaw. A hexagon key is included.

Code
994473

NuMedical



\$21.95
/box

25pcs/box


Plastic Bonding Trays

With adhesives to hold brackets in place. Can hold one full bracket set-up, with or without mixing wells.

Code	Colour
993430	D. Blue
993431	White
993432	Green
993433	Teal

Photography Mirrors


Rhodium-coated mirrors to assist dental professionals in taking intraoral photos.



Type	Price	Single-Sided
1	\$29.95 each	996131
2	\$29.95 each	996132
3	\$29.95 each	996133
4	\$29.95 each	996134
5	\$29.95 each	996135
6	\$29.95 each	996136
7	\$29.95 each	996137
8	\$29.95 each	996138

Photography Mirrors

Rhodium-coated mirrors to assist dental professionals in taking intraoral photos.



Type	Price	Single-Sided
9	\$29.95 each	996139
10	\$29.95 each	996124
11	\$29.95 each	996125
12	\$29.95 each	996126
13	\$29.95 each	996127
14	\$29.95/pack	996128

Type	Price	Double-Sided
15	\$39.95 each	996140
16	\$39.95 each	996141
17	\$39.95 each	996142
18	\$39.95 each	996143
19	\$39.95 each	996144

Orthodontic Retainer and Aligner Thermoforming Material, Hard

Made from a highly durable material that does not shrink crack or break. Offers excellent formability flexibility and superior strength. Ideal for fabrication of aligners and retainers. Also used to form other appliances such as splints nightguards positioners and temporary anterior bridges. Comes in square and round size.

Square size: 127mm (L) x 127mm (W).

Round size: 127mm (Diameter)

- BPA Free
- Non-Latex
- Does not bond to acrylics



SQUARE SIZE: 127mm (L) x 127mm (W)			
Thickness	Quantity	Price	Code
0.8mm	20pcs/bag	\$13.95/bag	994281
1.0mm	20pcs/bag	\$13.95/bag	994282
1.5mm	15pcs/bag	\$13.95/bag	994283
2.0mm	10pcs/bag	\$13.95/bag	994284



ROUND SIZE: 127mm (Diameter)			
Thickness	Quantity	Price	Code
0.8mm	20pcs/bag	\$13.95/bag	994285
1.0mm	20pcs/bag	\$13.95/bag	994286

Instructions:

1. Remove the pink protective sheet (Pic: 1). This side of the thermoforming polypropylene material should face upwards on the vacuum thermoforming machine.
2. Insert the thermoforming polypropylene material into the clamping frame (Pic: 2).
3. Place desired teeth model or cast on the vacuum plate (Pic: 3). It is recommended to spray the model with silicone spray for easy handling.
4. Allow 2-4 minutes for the thermoforming polypropylene material to heat sufficiently to achieve optimal sag for vacuum forming.
5. Thermoforming polypropylene material will be heated sufficiently when it achieves a sag between 15 - 20mm.
6. Once the material has been heated, lower the clamping frame over the model and activate the vacuum. It usually takes 10 - 20 seconds for thin materials, and 20 - 40 seconds for thicker materials to achieve the desired mold.
7. Switch off the vacuum thermoforming machine. Allow material to cool before handling. (Pic: 4)
8. Trim off excess material. (Pic: 5)
9. The mold is ready to be used.

Note: If bubbles start to appear on the thermoforming polypropylene material during the heating process, this is due to moisture. When this happens, peel off the protective sheets on both sides of the material. Store material in a sealed container with desiccants to fully remove the moisture. The material can be used after one week.

"BSA" Thermoforming/Coping Polypropylene Mouth Guard Material - SOFT

Made from a special formulated EVA plastic that is super soft and flexible. This material heats easily and conforms perfectly to the model. These custom-fit mouth guards allow for excellent oxygen intake.

*Note: To use, please remove the pink protective films.

- BPA Free
- Non-Latex
- Does not bond to acrylics





SQUARE SIZE: 127mm (L) x 127mm (W)

Thickness	Quantity	Price	Code
1.0mm	25pcs/bag	\$29.50/bag	994270
1.5mm	15pcs/bag	\$29.50/bag	994271
2.0mm	13pcs/bag	\$29.50/bag	994272
3.0mm	11pcs/bag	\$29.50/bag	994273



ROUND SIZE: 127mm (Diameter)

Thickness	Quantity	Price	Code
1.0mm	25pcs/bag	\$29.50/bag	994276
1.5mm	15pcs/bag	\$29.50/bag	994277
2.0mm	13pcs/bag	\$29.50/bag	994278
3.0mm	11pcs/bag	\$29.50/bag	994279



Materials from Germany, Hot in Japan and China.



Excellent Operation Experience




\$129.00
each

Distal End Cutter Standard

Hardness: HRC 60-65. Designed to cut the long distal end of wires up to 0.021" x 0.025" (0.53 mm x 0.64 mm). It securely holds the cut part.

Code
995911



\$129.00
each

Heavy Wire Cutter Standard

Hardness: HRC 60-65. Designed to cut various Australian wires, NiTi wires and stainless steel wires up to maximum diameter of 0.030" (0.76mm).

Code
995921



\$129.00
each

Distal End Cutter Long Handle

Hardness: HRC 60-65. Designed to cut the long distal end of wires up to 0.021" x 0.025" (0.53 mm x 0.64 mm). It securely holds the cut part.

Code
995912



\$129.00
each

Distal End Cutter Flush Standard

Hardness: HRC 60-65. Cuts and holds wire distally flush to the buccal tube or bracket while reducing debondings. Wires up to 0.021" x 0.025" (0.53 mm x 0.64 mm).

Code
995913



\$129.00
each

Light Wire Cutter 0° Standard

Hardness: HRC 60-65. Designed to cut ligature wires or elastic rings up to 0.021" (0.30 mm). Can't be used for cutting archwires.

Code
995909



\$129.00
each

Distal End Cutter Flush Long Handle

Hardness: HRC 60-65. Cuts and holds wire distally flush to the buccal tube or bracket while reducing debondings. Wires up to 0.021" x 0.025" (0.53 mm x 0.64 mm).

Code
995914

ORTHODONTICS

ORTHODONTICS

NuMedical

Code
995910

\$33.95
each



Light Wire Cutter 0° Long Handle

Hardness: HRC 60-65. Designed to cut ligature wires or elastic rings' up to 0.021" (0.30 mm). Can't be used for cutting archwires.

Code
995915

\$37.95
each



NiTi Cinch Back Plier Standard

Hardness: HRC 60-65. Designed to cinch back distal end of archwires to avoid causing injury to soft tissue. For archwires up to maximum diameter of 0.021" (0.53 mm).

Code
995916

\$54.95
each



NiTi Cinch Back Plier Long Handle

Hardness: HRC 60-65. Designed to cinch back distal end of archwires to avoid causing injury to soft tissue. For archwires up to maximum diameter of 0.021" (0.53 mm).

Code
995917

\$37.95
each




Light Wire Plier Standard

Hardness: HRC 60-65. Designed to bend archwires to make fine loops or spring loops.

Code
995918

\$33.95
each



Light Wire Plier Long Handle

Hardness: HRC 60-65. Designed to bend archwires to make fine loops or spring loops.

Code
995919

\$33.95
each



Light Wire Plier with Cutter Stand

Hardness: HRC 60-65. Designed to bend archwires to make fine loops or spring loops. The middle part of the jaw is used for cutting.

Code
995920

\$33.95
each



Light Wire Plier with Cutter Long Handle

Hardness: HRC 60-65. Designed to bend archwires to make fine loops or spring loops. The middle part of the jaw is used for cutting.

Code
995922

\$37.95
each



Heavy Wire Cutter Long Handle

Hardness: HRC 60-65. Designed to cut various Australian wires, NiTi wires and stainless steel wires up to maximum diameter of 0.030" (0.76mm).

Code
995923

\$33.95
each



V Stop Plier

Hardness: HRC 60-65. Designed to make V stop along archwires to prevent any movement.

ORTHODONTICS

ORTHODONTICS

NuMedical

Code
995924

\$33.95
each



Multi Function Plier

Hardness: HRC 52-55. Designed to bend, flatten and cut wires.

Code
995925

\$33.95
each



Torque Plier

Hardness: HRC 60-65. Designed to be used in pair to make torques on rectangle wires, thickness of tips: 0.050" (1.27mm).

Code
995926

\$33.95
each



Tweed Loop Plier

Hardness: HRC 60-65. Designed to make small loops, distal stops and Ω loops on archwires up to 0.022" (0.56mm).

Code
996753

\$43.95
each



Ortho Plier, Adam 64

Made with Tungsten Carbide inserts for cutting 0.7mm wires or less.

Code
996754

\$43.95
each



Ortho Plier, Adam 65 Flat/Round

Made with Tungsten Carbide inserts for cutting 0.7mm wires or less.

Code
996852

\$37.95
each



Triple Beak Pliers

Triple Beak Pliers

Code
996850

\$37.95
each



Tube Crimping Pliers

Tube Crimping Pliers

Code
996849

\$29.95
each



Reynolds Contouring Pliers

Reynolds Contouring Pliers 115

Code
996855

\$29.95
each



Ortho Bracket Applier

Ortho Bracket Applier

ORTHODONTICS

ORTHODONTICS

Nu²medical

Orthodontic Pliers

Code
995927

\$33.95
each



Band Remover

Hardness: HRC 52-55. Designed to remove bands and other bondable attachments in the posterior section.

Code
995928

\$136.75
each



Anterior Bracket Removing Plier

Hardness: HRC 52-55. Designed to remove anterior brackets.

Code
995929

\$136.75
each



Posterior Bracket Removing Plier Standard

Hardness: HRC 52-55. Designed to remove posterior brackets.

Code
995930

\$136.75
each



Posterior Bracket Removing Plier Long Handle

Hardness: HRC 52-55. Designed to remove posterior brackets.

Code
995931

\$37.95
each



Weingart Plier Standard

Hardness: HRC 52-55. Designed to cinch back distal end of archwires to prevent injury to soft tissue. NiTi wires' diameter up to 0.021" (0.53mm).

Code
995932

\$37.95
each



Weingart Plier Long Handle

Hardness: HRC 52-55. Designed to cinch back distal end of archwires to prevent injury to soft tissue. NiTi wires' diameter up to 0.021" (0.53mm).

Code
995933

\$37.95
each



Crimpable Hook Plier Standard

Hardness: HRC 52-55. Designed to fix the crimpable hooks along archwires.

Code
995934

\$136.75
each



Crimpable Hook Plier Long Handle

Hardness: HRC 52-55. Designed to fix the crimpable hooks along archwires.

Code
995935

\$136.75
each



Loop Bending Plier

Hardness: HRC 52-55. Designed to make different small loops for multi-loop edgewise archwires. Loop diameters: 2mm, 3.2mm and 4.2mm.

Orthodontic Pliers

\$136.75
each

Code
995936



How Plier

Hardness: HRC 52-55. Designed to hold archwires when guiding into buccal tubes, to hold ligature wires, or to form anterior band.


Code
995937



Kim Plier

Hardness: HRC 52-55. Designed to bend multi-loop edgewise wires, with cutter on tips.

Code
995938



Nance Loop Bending Plier

Hardness: HRC 52-55. Designed to make vertical loops of different height along archwires, 3.0, 4.0, 5.0, 6.0mm. Wires' diameter up to: 0.022" (0.56mm).


Code
995939



Bond Remover

Hardness: HRC 52-55. Designed to remove bond on the crown without any injury on enamel.

Code
995940



Three Jaw Plier

Hardness: HRC 52-55. Designed to make V loops and bend archwires. Wires' diameter up to 0.030" (0.76mm).


Code
995941



Omega Loop Bending Plier

Hardness: HRC 52-55. Designed to make fine Ω loops. Loops diameters: 0.063", 0.075", 0.10" (1.6mm, 2.0mm, 2.6mm).

Code
996146



Elastic Plier

Hardness: HRC 52-55. Designed to expand elastic separating rings before separating molar gap.

ORTHODONTICS

ORTHODONTICS



\$37.95
each

Bird Beak Pliers

Bird Beak Pliers

Code
996851



\$29.95
each

Crown Crimping Plier

To crimp the gingival margin of most sizes of aluminum and stainless steel crowns.

Code
996847



\$29.95
each

Crown Contouring Crimping Plier

To enhance crown form contours to improve gingival margins.

Code
996848

NuMedical

