

SURGICAL & KITS & IRRIGATION

SURGICAL & KITS & IRRIGATION

Nulimplant Drape Kit - 990726



	X2 Surgeon/Assistant Gowns	L
	X2 Surgical Cap	universal
	X4 Hand Towels	330mm x 230mm
	X1 Liquid Proof Mayo Cover	750mm x 1400mm
	X1 Absorbent/Liquid Proof Drape with Self-Adhesive Side	1000mm x 1500mm
	X2 Liquid Proof Drape	600mm x 600mm
	X1 Liquid Proof Drape	1200mm x 1200mm
	X1 Saliva Ejector	142mm x 6.5mm dia.
	X1 Surgical Aspirator Tips	Angled 30° 1/8"
	X1 Irrigation set for implant machine	universal
	X4 Sani Sleeve	1500mm x 55mm
	X6 Machine Protection Film	200mm x 200mm
	X4 Non-woven Sponge	75mm x 75mm
	X1 Biohazard Trash Bag	220mm x 355mm

\$183.00
/box

6 kits/box

Nulimplant Drape Kit

**First Implant Drape Kit Customized
for Australian Market.**

Code
990726

NuMedical

Nulimplant Kit 990706

1
BAG

\$35.95
/bag



	X3 Sani Sleeve	1500mm x 55mm
	X1 Liquid Proof Mayo Cover	750mm x 1400mm
	X2 Surgeon/Assistant Gowns	L
	X1 Absorbent/Liquid Proof Drape with Self-Adhesive Side	1000mm x 1500mm
	X1 Fluid-Repellent Fenestrated Drape with Adhesive O-Shape Cut and Loop	1350mm x 2000mm with 100mm dia.
	X1 Liquid Proof Drape	1600mm x 1900mm
	X4 Hand Towels	330mm x 450mm
	X1 Suction Tip	11mm dia.
	X1 Conical Adapter for Aspirators	0.6mm dia.
	X1 Saliva Ejector	142mm x 6.5mm dia.
	X2 Machine Protection Film	200mm x 200mm
	X1 Biohazard Trash Bag	220mm x 355mm
	X1 Transparent Bag for Physiodispenser or Control Unit	700mm x 600mm
	X1 Surgical Aspirator with Conical Connectors	2.2m

Nulimplant Kit

Code
990706

\$41.95
/bag

Nulplant Kit Plus 990717

1
BAG



Nulplant Kit Plus

	X2 Sani Sleeve	1500mm x 55mm
	X1 Liquid Proof Mayo Cover	750mm x 1400mm
	X2 Surgeon/Assistant Gowns	L
	X1 Surgical Drape	1000mm x 1000mm
	X1 Fluid-Repellent Fenestrated Drape with Adhesive O-Shape Cut and Loop	1350mm x 2000mm with 100mm dia.
	X1 Liquid Proof Drape	1600mm x 1900mm
	X2 Hand Towels	330mm x 450mm
	X1 Suction Tip	11mm dia.
	X1 Conical Adapter for Aspirators	0.6mm dia.
	X20 Sponge	75mm x 75mm
	X3 Machine Protection Film	200mm x 200mm
	X5 Cotton Balls	0.2g
	X1 Waste Collection Bag	910mm x 610mm
	X1 Surgical Aspirator with Conical Connectors	2.2m
	X1 Emesis Basin	250mm x 100mm x 50mm
	X4 Adhesive Tape	50mm x 200mm
	X1 Tube Holder	

Code
990717



Sterile Surgical Drape

ETO Gas Sterilised. Made of polypropylene with lamination. Ideal for dental implant surgeries to reduce bacterial transfer and contamination of the surgical site. These drapes are highly absorbent. Clean-up time is greatly reduced.

Size	Dimensions	Qty	Price	Code
1	450mm x 750mm	20pcs/box	\$13.95/box	990720
2	600mm x 600mm	20pcs/box	\$13.95/box	990721
3	900mm x 900mm	20pcs/box	\$19.95/box	990722



Sterile Surgical Drape with Hole

ETO Gas Sterilised. Ideal for dental implant surgeries to protect patients from dental polish fluids and blood during the surgical procedure. The triangular opening is designed for easy positioning and placement over the patient's face and mouth.

Size	Dimensions	Qty	Price	Code
2	600mm x 800mm	10pcs/box	\$13.95/box	990710
3	900mm x 1200mm	10pcs/box	\$19.95/box	990711

Nu
medical

SURGICAL & KITS & IRRIGATION

SURGICAL & KITS & IRRIGATION



\$18.95
/box

3pcs x 120packs/box

Sterile Non-Woven Gauze

Made from 70% rayon and 30% polyester which absorption capacity is more stronger than regular cotton gauze materials. Steril package saves time. Each pack comes with 3 pieces and 120 packs in box. EO Gas sterile. Disposable.
Size: 75mm x 75mm.

Code
992817



\$18.95
/box

5pcs x 120packs/box

Sterile Cotton Rolls

Made from 100% cotton which absorption capacity is more stronger. Steril package saves time. Each pack comes with 5 pieces and 120 packs in box. EO Gas sterile. Disposable.
Size: 10mm x 38mm.

Code
992780



Sterile Cotton Tip Applicators

Made from 100% cotton which absorption capacity is more stronger. Steril package saves time. Each pack comes with 2 pieces and 100 packs in box. EO Gas sterile. Disposable.

Size	Qty	Price	Code
80mm (L)	2pcs x 100packs/	\$8.90/box	993766
150mm (L)	2pcs x 100packs/	\$9.90/box	993767



\$99.95
/box

50pcs/box

Sterile Pressure Pack with TIE

EO Gas sterile. Single-use. 2 pieces per steril package and 50 packages in a box.

Code
990714



\$16.50
/box

50pcs/box

Sterile Barrier Film

Disposable recyclable sleeves for handpieces and hoses. Sterilised and packed individually. Single-Use only. Size: 1500mm (L) x 55mm (W) x 30mm (Diameter)

Code
990723



\$6.95
/bag

knitted cuffs



sewn by ultra-sonic technology

Sterile Gown Kit

Made of heavy SMS material, breathable, comfortable and reliable protection apparel with knitted cuffs. Each kit contains 1 x Surgeon Gown and 2 x Hand Towel. Seams are sewn by ultra-sonic technology. Latex-free and fluid-resistant. Size: L 1500mm x W 1200mm

Code
991538

SURGICAL & KITS & IRRIGATION



\$29.95
/box

Surgical Aspirator Tips

Made from 100% durable polyurethane latex-free. Fits 11mm suction holders. EO Gas sterile. Single-use.

Code
990703



\$41.50
/box

50pcs/box

Sterile Aspirator Tips

With ergonomic hand-grip and aspiration regulator. EO Gas sterile. Single-use.

Code
990709



\$49.95
/box

10pcs/box

Sterile Irrigation Set

Ideal for surgical units. EO Gas sterile. Single-use. 990716 is ideal for brand implant machine, like NSK/WH/Saeyang/Saeshin/Nobelbircare/Dentium/Staumaun.

Code Size
990716 Universal



\$49.95
/box

10pcs/box

Sterile Irrigation Set


Ideal for surgical units. EO Gas sterile. Single-use. 990737 is ideal for brand implant machine of Kavo.

Code Size
990737 Kavo



\$27.25
/box

10pcs/box



Sani Sleeves

Disposable recyclable sleeves for handpieces and hoses. Sterilised and packed individually. Single-Use only. Size: 1500mm (L) x 55mm (W) x 30mm (Diameter)

Code
992491



\$35.95
/box

8 kits/box

Sterile Suction Tube with Connectors

Sterile package for saving time. 2.2m length. EO Gas sterile. Single-Use.

Code
990707