

Titanium Scaler Tips

Titanium Scaler Tips



\$37.95
/pack

5pcs/pack

NuMedical Scaler Tips GD1T

Supragingival scaling. Compatible with SATELEC and DTE.

Code
995827

Surgery Tips, Stainless Steel

NuMedical

ENDODONTICS TREATMENT			
TYPE	Price	Code	
UE1, Diamond-Coated For apical root debridement *Diamond grit 40um *MECTRON EN1*	\$72.10	995794	
UE2 For gentle apical root debridement *MECTRON EN2*	\$72.10	995795	
UE3, Diamond-Coated For apical root debridement *Diamond grit 40um *MECTRON EN3*	\$72.10	995796	
UE4 For gentle apical root debridement *MECTRON EN4*	\$72.10	995797	

Titanium Scaler Tips

Surgery Tips, Stainless Steel



BONE SURGERY			
TYPE	Price	Code	
US1 For osteotomy *MECTRON OT6*	\$105.95	995768	
US1L For osteotomy *MECTRON OT8L*	\$105.95	995769	
US1R For osteotomy *MECTRON OT8R*	\$105.95	995770	
US2 For osteotomy *MECTRON OT7*	\$105.95	995771	
US3 For precise osteotomy in thin/narrow regions *MECTRON OT2*	\$105.95	995772	
US4 For bone chip harvesting and crown lengthening *MECTRON OP3*	\$105.95	995773	
US5 For bone chip harvesting and bone modelling *MECTRON OP1*	\$105.95	995774	
US6 For precise osteotomy *MECTRON OT7S-3*	\$105.95	995775	
SINUS LIFT SURGERY			
	Price	Code	
UL1 For peeling of the sinus membrane *MECTRON EL3*	\$105.95	995776	
UL2 For peeling of the sinus membrane *MECTRON EL1*	\$105.95	995777	
UL3, Diamond-Coated Non-traumatic bone cutting on thin bone region *MECTRON OT5*	\$105.95	995778	
UL4 For micrometric osteotomy near soft tissue *MECTRON OT1*	\$105.95	995779	
UL5 95° angle, for peeling of the sinus membrane	\$105.95	995780	
EXCELCYMOSIS			
	Price	Code	
UC1 For cutting ligaments, with minimal wound *MECTRON EX1, PIEZO LC1*	\$105.95	995786	







Titanium Scaler Tips

NuMedical

Titanium Scaler Tips

Titanium Scaler Tips

IMPLANTS			
	Price	Code	
UI1 For implants with tip ball Ø1.6mm	\$105.95	995781	  
UI2 For implants with tip ball Ø2.0mm irrigation from tip point	\$105.95	995782	  
UI7, Diamond-Coated To finalise the implant site and for bone preparation *Diamond grit: 100um *MECTRON OT4*	\$105.95	995783	  
UI8 For implants, extra long with tip ball Ø1.6mm	\$105.95	995784	  
UI9 For implants with tip ball Ø2.8mm	\$105.95	995785	  

PERIODONTAL SURGERY			
	Price	Code	
UP1 For root scaling *MECTRON PS1*	\$97.25	995787	  
UP2 For root scaling *MECTRON PS2*	\$97.25	995788	  
UP3 For root scaling *MECTRON PS6*	\$97.25	995789	  
UP4, Diamond-Coated For root surface debridement and root planning *Diamond grit 40um *MECTRON OP5*	\$97.25	995790	  
UP5 For perio treatment *MECTRON PP1*	\$97.25	995791	  
UP6, Diamond-Coated For root preparation in perio surgery *Diamond grit 40um *MECTRON OP6*	\$97.25	995792	  
UP7, Diamond-Coated For bone remodeling and root canal enlargement *Diamond grit 40um	\$97.25	995793	  


Titanium Scaler Tips


Scaler Tips, Stainless Steel (TITANIUM COATING)


SCALING			
TYPE	Price	Code	FOR WOODPECKER & EMS SCALER
G1T Supragingival Scaling	\$13.75	995712	  
G2T Supragingival Scaling <i>*EMS B*</i>	\$13.75	995714	  
G3T Supragingival Scaling <i>*EMS P*</i>	\$19.00	995716	  
G4T Supragingival Scaling <i>*EMS A*</i>	\$19.00	995718	  
G5T Supragingival Scaling <i>*SATELEC 1, MECTRON S1, NSK G1*</i>	\$19.00	995720	  
G6T Supragingival Scaling <i>*EMS C, SATELEC 2, MECTRON S3, NSK G2*</i>	\$19.00	995722	  
G7T Crown Removal <i>*EMS D, SATELEC 5, NSK G25*</i>	\$25.80	995724	  
G8T Crown Removal with irrigation <i>*SATELEC 5AE, NSK G26*</i>	\$25.80	995726	  

SCALING			
TYPE	Price	Code	FOR DTE & SATELEC SCALER
GD1T Supragingival Scaling	\$13.75	995713	  
GD2T Supragingival Scaling <i>*EMS B*</i>	\$13.75	995715	  
GD3T Supragingival Scaling <i>*EMS P*</i>	\$19.00	995717	  
GD4T Supragingival Scaling <i>*EMS A*</i>	\$19.00	995719	  
GD5T Supragingival Scaling <i>*SATELEC 1, MECTRON S1, NSK G1*</i>	\$19.00	995721	  
GD6T Supragingival Scaling <i>*EMS C, SATELEC 2, MECTRON S3, NSK G2*</i>	\$19.00	995723	  
GD7T Crown Removal <i>*EMS D, SATELEC 5, NSK G25*</i>	\$25.80	995725	  
GD8T Crown Removal with irrigation <i>*SATELEC 5AE, NSK G26*</i>	\$25.80	995727	  

Titanium Scaler Tips

SCALING			
TYPE	Price	Code	FOR AMDENT SCALER
A1T Supragingival Scaling *AMDENT 37*	\$35.60	995710	  
A2T Supragingival Scaling *AMDENT 39*	\$35.60	995711	  




PERIODONTAL			
TYPE	Price	Code	FOR WOODPECKER & EMS SCALER
P1T Periodontal Scaling *MECTRON S2*	\$25.80	995728	  
P3T Periodontal Scaling *EMS PS, MECTRON P10*	\$25.80	995730	  
P4T Periodontal Scaling *NSK G4*	\$25.80	995732	  
P2LT Periodontal Scaling *NSK P25L, AMDENT PE-40L*	\$38.50	995734	  
P2RT Periodontal Scaling *AMDENT PE-40R*	\$38.50	995736	  


PERIODONTAL			
TYPE	Price	Code	FOR DTE & SATELEC SCALER
PD1T Periodontal Scaling *MECTRON S2*	\$25.80	995729	  
PD3T Periodontal Scaling *EMS PS, MECTRON P10*	\$25.80	995731	  
PD4T Periodontal Scaling *NSK G4*	\$25.80	995733	  
PD2LT Periodontal Scaling *NSK P25L, AMDENT PE-40L*	\$38.50	995735	  
PD2RT Periodontal Scaling *AMDENT PE-40R*	\$38.50	995737	  

Titanium Scaler Tips

Titanium Scaler Tips

NuMedical

ENDODONTICS			
TYPE	Price	Code	FOR WOODPECKER & EMS SCALER
E1T 120° needle holder, for root canal cleaning *NSK E11*	\$29.65	995738	
E2T 95° needle holder, for root canal cleaning *NSK E12*	\$29.65	995740	
E3T For lateral thermal condensation of G/Percha *EMS H, SATELEC S07, MECTRON D3, NSK E6*	\$29.65	995742	
E4T Removes obstructions in the root canal *SATELEC ET40*	\$42.35	995744	
E5T Removes fillings in the coronal third *SATELEC ET20*	\$42.35	995746	
E8T Bur Holder, used to expand root canal and grind teeth	\$34.15	995748	
E9T Bur Holder, used to expand root canal and grind teeth *NSK V30*	\$34.15	995750	
E10T For root canal treatment *MECTRON EN1*	\$34.15	995752	
E11T For root canal treatment *MECTRON EN3*	\$34.15	995754	

ENDODONTICS			
TYPE	Price	Code	FOR DTE & SATELEC SCALER
ED1T 120° needle holder, for root canal cleaning *NSK E11*	\$29.65	995739	
ED2T 95° needle holder, for root canal cleaning *NSK E12*	\$29.65	995741	
ED3T For lateral thermal condensation of G/Percha *EMS H, SATELEC S07, MECTRON D3, NSK E6*	\$29.65	995743	
ED4T Removes obstructions in the root canal *SATELEC ET40*	\$42.35	995745	
ED5T Removes fillings in the coronal third *SATELEC ET20*	\$42.35	995747	
ED8T Bur Holder, used to expand root canal and grind teeth	\$34.15	995749	
ED9T Bur Holder, used to expand root canal and grind teeth *NSK V30*	\$34.15	995751	
ED10T For root canal treatment *MECTRON EN1*	\$34.15	995753	
ED11T For root canal treatment *MECTRON EN3*	\$34.15	995755	